



# The *mar345* Software Guide

Version 1.2  
April 4, 2001

Written by Dr. Claudio Klein

Copyright 1997-2001 X-Ray Research G.m.b.H.

**All Rights Reserved.**

This document contains proprietary and confidential information of X-ray Research G.m.b.H. The contents of this document may not be disclosed to third parties, copied or

## 1. Distribution Media

Media	OS	Suggested directory	Mount CD-Rom with
CD-Rom	IRIX 5.x,6.x	/usr/people/mar345	mount -t iso9660 -r /dev/rdisk/dksXdYvol /CDROM
	DEC Unix 5.0	/usr/users/mar345	mount -t iso9660 -r -o rrip /dev/rzYc /mnt
	Linux, 2.0.x	/home/mar345	mount -t iso9660 -r /dev/cdrom /mnt/cdrom

X denotes the SCSI-controller (e.g.0) and Y the SCSI-id (e.g. 4).

Most executable programs are also available as gzip compressed files via anonymous FTP from ftp.EMBL-Hamburg.de, directory pub/marx/mar345/.... Due to their size (180 MB), scanner calibration files are not available per FTP. Please contact X-ray Research for assistance.

## 2. Compilation Notes

OS	Compilers used	Motif version used	
IRIX	C-7.1	Motif 1.2.4	X11R6

## 3. Environment

The following logicals assignments must be set to run certain programs:

Variable	Description	Used by
MARHOME	Master directory of software distribution.	-
MARLOGDIR	Directory where log files will be created.	mar345, mar345xf, scan345, marHKL

## 4. Directory Structure

The *mar345* software distribution directory (\$MARHOME) contains the following subdirectories:

Subdirectory	Contents
bin/sgi bin/osf bin/linux	Binary executables for corresponding platforms.
man/1 man/cat1 man/man1	Unformatted man pages for selected programs. Formatted man pages (SGI). Compressed unformatted man pages (Linux & DEC/Linux)
man/ps man/doc man/help	Postscript files of formatted man pages, ready for printing. ASCII text of formatted man pages, ready for online read (more). Online help files for GUI's ( <i>mar345</i> , <i>marView</i> , <i>marTools</i> ).
log log/log log/lp log/spy log/sets	Log-files for programs <i>mar345</i> and <i>scan345</i> . Up to 99 versions of mar.log files. Up to 99 versions of mar.lp files (statistical output). Up to 99 versions of mar.spy files (native controller messages). Data collection strategy template files for <i>mar345</i> .
src tables	Source code of selected programs. Scanner specific calibration and configuration files.
Optional: marfilm	Latest <i>marFLM</i> distribution.
marxds	Latest <i>marXDS</i> distribution.

## 5. Description of Programs

Name	man page	Priority	Description
<b>mar345</b>	mar345	A	GUI for data collection and display.
<b>marstart</b>	-	A	Works together with program <i>mar345</i> and must be available.
<b>marView</b>	marView	A	GUI for data display.
<b>marTools</b>	marTools	A	GUI for image file manipulation: format conversion, etc.

<b>marHKL</b>	-	B	(not installed in a standard SGI installation). GUI for HKL processing package (DENZO, xdisp).
<b>marFLM</b>	-	B	GUI for MOSFLM processing package.
<b>marXDS</b>	-	B	GUI for XDS processing package (XDS, xscale).

Hardware related programs:

The documentation can be found in directory \$MARHOME/man. Several formats are available:

Directory	Description
man/1	Unformatted man pages for selected programs.
man/cat1	Formatted man pages (SGI).
man/man1	Compressed unformatted man pages (Linux & DEC Unix).
man/ps	Postscript files of formatted man pages, ready for printing.
man/doc	ASCII text of formatted man pages, ready for online read ( <i>more</i> ).
man/help	Online help files for GUI's ( <i>mar345, marView, marTools</i> ).

Name	Description
<u>Programs:</u>	
mar345	Documentation for program <i>mar345</i> .
marView	Documentation for program <i>marView</i> .
marTools	Documentation for program <i>marTools</i> .
tkmarcvt	Documentation for program <i>tkmarcvt</i> .
mar345xf	Documentation for program <i>mar345xf</i> .
scan345	Documentation for program <i>scan345</i> .
marsim	Documentation for program <i>marsim</i> .
marcvt	Documentation for program <i>marcvt</i> .
spiralpack	Documentation for program <i>spiralpack</i> .
catmar	Documentation for program <i>catmar</i> .
marPeaks	Documentation for program <i>marPeaks</i> .
marIndex	Documentation for program <i>marIndex</i> .
marPredict	Documentation for program <i>marPredict</i> .
marStrategy	Documentation for program <i>marStrategy</i> .

Others:

mar345\_formats Documentation for *mar345* image formats.

## 7. Installation on SGI IRIX

### 7.1 Installing the Software (IRIX 6.5)

#### 7.1.1 Login as user "root" and create a new user account called **mar345**

Do the following:

- a) From the toolchest (usually in the upper left corner of the screen) choose "**System Manager**" from the "**System**" menu.
- b) In the main window of the "**System Manager**", choose: "**Security and Access Control...**"
- c) Select "**Add a User Account**" and go through steps 1 to 10.
  - Choose "**mar345**" as "**Login name**" (suggested).
  - Choose "**Password**", "**User-ID**" and "**Primary Group**" as desired.
  - Choose "**/usr/people/mar345**" as "**Home Directory**" (suggested).
  - Choose "**tcsh**" as "**Shell Program**" (strongly recommended).
  - As a last step, press "**OK**" to accept the settings.

#### 7.1.2 Login as user: **mar345**

#### 7.1.3 Installing the software from CD-ROM:

Insert the CD-ROM in the CD-ROM reader. After a couple of seconds, the automounter should mount the CD automatically as **/CDROM**. If not, the super user should try:

```
mount -t iso9660 -r /dev/rdisk/dks0d4vol /CDROM
```

where 4 is the SCSI-unit number of the CD-ROM on the system, which may vary. If you are not certain, which SCSI-unit your CD-ROM drive uses, try command:

```
hinv
```

to get a "hardware inventory" of your system which will say something like:

```
CDROM: unit 4 on SCSI controller 0:
```

When successful, the contents of the CD-ROM should be copied into the login directory of the account **mar345**. To do so, as user "**mar345**" type:

```
/CDROM/mar_install
```

The installation script chooses reasonable defaults that may be accepted or modified.

## 7.1.2 Setting Up the Ethernet Connection

The scanner control program *mar345* communicates with the scanner through a Ethernet interface. To use *mar345*, the network must be configured to meet the requirements of the mar controller. The mar scanner will respond only to requests made to IP-address 192.0.2.1. The mar scanner also requires the host computer Ethernet interface to be set to address 192.0.2.2. These addresses belong to a pool of addresses that are not assigned to official networks so there should not be any conflict with the outside world. Do the following:

7.1.2.1 Edit file `/etc/hosts`, e.g. type:

```
jot /etc/hosts
```

```
192.0.2.1    mar345
192.0.2.2    gate-$HOSTNAME
192.0.2.2    gate-mars
```

where `$HOSTNAME` is the original name of your computer, e.g. if the computer is called "mars", then write:

7.1.2.2 If the second network interface is not yet configured, edit file `/etc/config/netif.options`:

```
jot /etc/config/netif.options
```

and modify lines:

```
: if2name=
: if2addr=gate-$HOSTNAME
```

to:

```
if2name=ec1
if2addr=gate-$HOSTNAME
```

7.1.2.3 Reboot.

7.1.2.4 Check settings

The second network interface (`ec1`) should now be configured and you should be able to access other hosts (e.g *mar345*) on network 192.0.2. To check, type:

```
/usr/etc/ifconfig ec1
```

This command should come back with something like:

## 8. Installation on DEC Unix

### 8.1 Installing the Software (DEC Unix 4.0)

8.1.1 Create a new user account called *mar345*.

Type:

a) **su**

= become superuser

b) `type /bin/ADM4 (then press enter)`

= "The "Account Manager" graphical interface

c) In the main window of the "Account Manager", choose: "Add..."

d) In the window "Create/Modify Local User" fill in at least the "User name" and "Shell" press "OK".

As default shell, the "/bin/tcsh" is recommended. However, the "tcsh" is not installed by default on DEC Unix 4.0. If "tcsh" is unavailable (look also in /usr/bin and /usr/local/bin) choose: "/bin/csh".

As "Home Directory" use of **/usr/users/mar345** is suggested. At this time you may create a

... for this user account but you don't have to. All other fields do not need to be modified

## 8.2 Setting Up the Ethernet Connection

The scanner control program *mar345* communicates with the scanner through a Ethernet interface. To use *mar345*, the network must be configured to meet the requirements of the mar controller. The mar scanner will respond only to requests made to IP-address 192.0.2.1. The mar scanner also requires the host computer Ethernet interface to be set to address 192.0.2.2. These addresses belong to a pool of addresses that are not assigned to official networks so there should not be any conflict with the outside world.

Do the following:

### 8.2.1 Configure the second Ethernet card (tu1)

Type:

- 1.) **su** = become superuser
- 2.) **/usr/sbin/netconfig** = call the "Network Configuration" GUI
- 3.) Choose **"Ethernet Interface – tu1"**
- 4.) **Select: "tu1 – Interface Configuration:"** **Enable**
- 5.) Enter **"Host name:"** **marcontrol**
- 6.) Enter **"IP address:"** **192.0.2.2**

```
vi /etc/hosts
```

and add the following line to the end of the file:

```
192.0.2.1    mar345
```

If you can't find an entry for name marcontrol, also add:

```
192.0.2.2    marcontrol
```

## 9.1 Installing the Software

### 9.1.1 Create a new user account called *mar345*.

Type:

- a) **su** = become superuser
- b) **adduser mar345** = adds user account "mar345" using defaults.
- c) Edit **/etc/passwd** and change the default shell of user "**mar345**" to **/bin/tcsh**.  
Possibly, the password of account "**mar345**" has to be set.

### 9.1.2 Login as user: **mar345**

### 9.1.3 Installing the software from CD-ROM:

1. Insert the CD-ROM into the CD-ROM reader. When the system prompts the CD-ROM is going to be

**9.2.1** Configure the second Ethernet card (eth1)

- Type:
- 1.) **su** = become superuser
  - 2.) **/usr/bin/netcfg** = call the "Network Configuration" GUI
  - 3.) Choose section "**Interfaces**"
  - 4.) Choose button "**Add**"
  - 5.) Choose "**Interface Type:**" **Ethernet**
  - 6.) In window "**Edit Ethernet/Bus Interfaces**"
    - Enter: "**IP**": **192.0.2.2**
    - Enter: "**Netmask**": **255.255.255.0**
    - Select: "**Activate interface at boot time**"
  - 8.) Press: **Done**
  - 9.) Confirm: **Save configuration**
  - 10.) In main window, select interface "**eth1**" and click button "**Activate**"
  - 11.) Click button "**Save**" and "**Quit**".

For the kernel to actually initialize 2 Ethernet cards, please consult the documentation of the Linux boot load "**lilo**", i.e. read the lilo man page or the lilo documentation in `/usr/doc/lilo*`.

**Note:** On SuSE systems, the integration of the second Ethernet card is done via the SuSE system administration tool (`yast` or `yast2`). Program `netcfg` is not available.

**9.2.2** Edit file `/etc/hosts`, e.g. type:

```
vi /etc/hosts
```

and add the following line to the end of the file:

```
192.0.2.1    mar345
```

If you can't find an entry for name `marcontrol`, also add:

```
192.0.2.2    marcontrol
```

The second network interface (eth1) should be configured now and you should be able to access other hosts (e.g. `mar345`) on network `192.0.2`. To check, type:

```
ifconfig eth1
```

This command should come back with something like:

```
eth1      Link encap:10Mbps Ethernet  HWaddr 00:80:C6:FF:EF:08
          inet addr:192.0.2.2  Bcast:192.0.2.255  Mask:255.255.255.0
          UP BROADCAST RUNNING MULTICAST  MTU:1500  Metric:1
          RX packets:0 errors:0 dropped:0 overruns:0
          TX packets:0 errors:0 dropped:0 overruns:0
          Interrupt:12 Base address:0x320
```

The correct routing table can be checked using command:

```
netstat -r
```

It should say something like:

**Kernel IP routing table**

Destination	Gateway	Genmask	Flags	MSS	Window	irrt	Iface
193.141.161.0	*	255.255.255.0	U	1500	0	0	eth0
192.0.2.0	*	255.255.255.0	U	1500	0	0	eth1
127.0.0.0	*	255.0.0.0	U	3584	0		

**9.2.3** Connect the scanner to the second Ethernet card and power it up.

## 10. Running mar345

10.1 Turn the scanner on. Wait until you hear the plate locking sound (click!).

10.2 Login as user mar345.

10.3 Try to talk to the scanner, i.e. type:

```
ping mar345 or /usr/etc/ping mar345 (SGI)
```

If the scanner is accessible, ping comes back with:

```
PING mar345 (192.0.2.1): 56 data bytes  
64 bytes from 192.0.2.1: icmp_seq=0 ttl=255 time=1 ms  
...
```

If ping comes back with:

```
ping: mar345: Unknown host
```

then, mar345 has not been inserted into file /etc/hosts (see 2.1).

If ping hangs with:

```
PING mar345 (192.0.2.1): 56 data bytes
```

then the reason might be  
a) the interface has not been configured correctly.