

# Pixirad-1: Area Detector with CdTe Sensor Detector Pool Guide







## **Quick Start: Mechanical Connections**

Air flow, OD: 4mm Necessary to avoid condensation within the detector enclosure. Flow rate: 0.5 liters/minute ((≤3 liter/min for best performance).

NOTE: Please disconnect the air from the detector while adjusting the flow. Over pressuring will tear the Mylar film!



Cat 5 or 6: Connects to 2<sup>nd</sup> ethernet card on PC (in expansion bus)

# **Quick start**

- Mounting detector to your experimental setup:
  - Please mount with M3 screws several mounting holes are available on the sides or bottom (remove rubber feet) of the detector head.
  - Mounting plate is available from the DP staff.
- IOC Startup:
  - Login information:
    - Ask Detector Pool for Login information. Ext 2-9490
  - Ensure that the chiller is connected and running
  - Ensure that the detector power supply is turned on (switch on back of 2<sup>nd</sup> box).
  - IOC startup icon is on the desktop:
    - Clicking will launch both the ioc and medm screen



- After dry gas has flowed for a few minutes, turn cooling on. Recommended to run at -30°, but ensure that setting is not below dew point or sensor will ice up.
- HV settings: Run at 400V once chilled. Leave "HV Mode" in "auto" and "HV State" in "off". HV is turned on automatically once you start an acquisition.
- Run "Auto Calibrate" before taking data.

3

#### **Medm Screen**

🔀 pixirad.adl	B		
Pixirad Detector Control - 13PR1:cam1:			
Setup	Shutter	High Voltage	
asyn port PIXI	Shutter mode None 🖃	HV mode Auto 🖃 Auto	
EPICS name 13PR1:cam1:	Status: Det. Closed EPICS Closed	HV state <mark>_0ff </mark> ⊒]Off	
Manufacturer Pixirad	Open/CloseOpenClose _	HV setpoint 400.0 V	
Model Pixirad 1	Delay: Open 0.000 Close 0.000	HV value 0.0V	
Connected	EPICS shutter setup 📃 🖳	HV current 0.0 nA	
Connection Connect Disconnect	Collect	Environmental	
Debugging 🖳		Cooling 💁 🖉 Off	
Plugins	Acquire period 0.000	Setpoint 5.0C	
	# Tmages 1000 1000	Cold Temp. 22.8C	
Stats 9 Other 9	Images collected 0	Hot Temp. 17.5C	
	Colors collected 0	Box Temp. 24.8C	
Detector	UDP buffs read()	Box humidity 14.5%	
$\begin{array}{c} \text{Detector Size } 470 & 512 \\ \text{Threehold 1 (keV) } 10.0 & 10.0 \\ \end{array}$	UDP speed (MB/s)0.0	Dew point -4.1C	
Threshold 2 (keV) $\frac{10.0}{20.0}$ 29.7	UDP buffs free/max 1500/1500	Peltier power 0.0%	
Threshold 3 (keV) $50.0$ 29.7	Frame type 2 color DTF 2 Color DTF	Status UK	
Threshold 4 (keV) $30.0$ 29.7	Trigger mode Internal Internal	Attributes	
Sunc in polarity Pos,	Done	File	
Sync out polarity Post POSt	Acquire Start Stop		
Sync out function Read done - Read dor	Image counter p		
Auto calibrate AutoCal	Array callbacks Frable of Frable		
System reset Reset Done			
	Status		
	Status: Server returned OK		
	To server: DAQ:! INIT 5.0 0 400.0 0		
From server: DETECTOR 1022 GDT: DAQ:! INIT 5.0 0 400.0 0%r			

## A few notes about detector operation

- Detector has 2 counters per pixel, and two discriminators per counter
  - Thresholds determine discriminator settings
  - Frame type field configures discriminators and counters:
    - 1 color low acquires single image using lower discriminator
    - 1 color high acquires single image using upper discriminator (threshold 2)
    - 2 color implements both discriminators in single counter. Actually records 2 images to output file.
    - DTF modes Reads out single counter while other counter is being used to acquire data. Permits fastest data acquisition... Not well tested!
      - 2 color DTF uses threshold 1 and 3
- Images are saved through Area Detector plug-ins. 2 or 4 color modes can only be saved through NetCDF or HDF formats.
- First image of multi-image acquisitions is always blank.
- "Stop" image acquisition does not work properly, and will require detector reset (red button on lower left) before resuming detector operations.
- "Acquisition Period" is used to increase the time between images. There is a 8ms readout time per counter, and the detector will wait before starting the next acquisition. To maximize your framerate for a given exposure time, set "Acquisition Period = 0.0"

#### Viewers

- An ImageJ shortcut is available on the desktop
- If ImageJ fails to display your images:
  - Double check that you have the detector properly named in the ImageJ plugin.
  - On the medm screen, ensure that "array callbacks" and the Image1 plug-in are both enabled.





#### The Chromatic Photon Counting

Ó.

Web page and content by: Gloria Spandre

#### The PiXirad imaging sensors

PIXIRAD-1: a single unit system, 250K pixels, 500K counters 3×2.5 cm<sup>2</sup> active area

Sensor specs:	CdTe, 650 µm, 30.9 × 25.0 mm² Schottky type diode Electron collection at pixel
ASIC+CdTe base block	512 × 476 pixels
Number of blocks	1
Global active area	31 x 25 mm <sup>2</sup>
Total number of pixels	243712
Total number of counters	487424
Pixel size	60 µm hexagonal arrangement
Pixel density	323 pixels/mm <sup>2</sup> , equivalent to 55 µm on square arrangement
Pixel rate capability	10 <sup>6</sup> counts/pixel/s (after dead-time correction)
Global rate capability	2.4x10" counts/s
Pixel dead-time	300 ns
Position resolution	11 line pairs/mm at MTF 50%
Reading while taking data	possible
Energy range	1-100 keV
Detection efficiency @10 keV, 25 keV, 50 keV	100%, 100%, 98%
Counters depth	15 bits
Read-out time @50 MHz clock	5 ms/counter
Frame rate	200 readouts/s
Minimum applicable global threshold	200 electrons
Sensor bias voltage	200 + 400 V
Leakage current density	5 nA /cm <sup>2</sup> at 400 V, -20 °C
Typical number of defective pixels	less than 1%
Number of independent thresholds (colors)	2 set of two (swappable in real time)
Camera specs:	
Size (W×L×H)	14×14×7 cm <sup>3</sup>
Weight	< 2Kg
Power consumption	60 Watts (typical)
Cooling	liquid or forced air
Operating temperature	+40 -40 °C