

Vortex 4-element SDD Quick Start Guide

1. **Connect fiber cable** from PC to mini-PXI crate
2. **Turn on the mini-crate** (PXI/CompactPCI) with XMAP module.
3. Turn on (or reboot) & **log-in** as “This computer”:
 - User: *lang* (or *dxp*)
 - Password: ******* (*Ask DP Staff*)
4. **Start EPICS IOC**
 - Double-Click on desktop icon: “*dxp_iocXmapDP1*”
 - EPICS PV prefix: **dxpXmapDP1**:
5. **Start MEDM** control screens
 - Double-click on desktop icon: “*dxp_MEDM*”
6. **MEDM screens from APSshare**
 - You can also access the medm screens from your local beamline computer if you mount APSshare from the dserv’s.
 - The medm startup scripts are located in:
/APSshare/DetectorPool/ For example:
 - *start_vortex_4element_xmap_dxpXmapDP1*
 - *start_vortex_4element_xmap_dxpXmapDP2*
 - *start_vortex_4element_xmap_dxpXmapDP3*



WARNING: Ion Pump AC adapter or power supply must be on continuously!!