


Using an Amptek (Fluorescence-) Detector at Beamline 27-ID

05/17/2016, tg

1. Hardware Setup



Connect


- power cable **Detector** ↔ **Detector Controller** “POWER”
- signal cable **Detector** “OUT” ↔ **Detector Controller** “AMP IN”
- signal cable(s) **Detector Controller** “AMP OUT” ↔
 - **Canberra ADC** “ADC IN” (for full MCA spectrum / “soft” counting within regions of interest)
 - **Ortec SCA** “IN” (for “hard” counting within a region of interest)
 - **Oscilloscope** (for trouble shooting / calibration)

Note: Combination of devices connected to **Detector Controller** “AMP OUT” (ADC, SCA(s), Oscilloscope) can be freely chosen, but will affect calibration of the detector. Therefore, “AMP OUT”-configuration should not be changed after the detector is calibrated.

2. Acquiring Full MCA Spectra


2.1 EPICS Setup

- Canberra ADC setup: Open MEDM window for ADC



- Configure ADC as shown above
- Use "Gain/Range" setting for total number of channels and resolution (#channels/10 V)
- Use "Lower (Upper) Level Disc." to suppress noise at lower (upper) end of spectrum

- **MCA setup:** Open MEDM window for MCA and MCA Setup




MCA Setup	
Max. number of channels	2048
Number of channels to use	1450
Number of channels read	1450
Preset Live Time	0.00
Preset Real Time	0.00
Preset Start Channel	0
Preset End Channel	0
Preset Counts	0
Dead Time Warning (%)	40.00
Dead Time Error (%)	70.00
Mode	PHA
Channel Advance	Internal
MCS Int. Dwell Time	0.00e+00
MCS Ext. Prescale	1
Read spectrum every	1 second
Check done every	.1 second
Wait for client	Disable
Client wait	Done

- set number of channels to include
- set "Live" or "Real" acquisition time, 0.00 for continuous acquisition

Plot can be customized through plot setup menu (right-click in plot window)

Regions of interest (ROIs) can be defined in the extended MCA window, accessible through the “More”-menu



X Axis

Axis Style: linear

Axis Range: user-specified

Minimum Value: 0.000000

Maximum Value: 1600.000000

Time format: hh:mm:ss

Y1 Axis

Axis Style: linear

Axis Range: auto-scale

Minimum Value: .

Maximum Value: .

Y2 Axis

Axis Style: linear

Axis Range: auto-scale

Minimum Value: .

Maximum Value: .

Close

Note: For the Canberra ADC/AIM modules to be recognized by EPICS, these modules have to be fully booted-up before the associated IOC “27idgen1” is booted. If Canberra modules are not recognized, then **(STAFF only!)**

- Turn-off IOC “27idgen1” VME crate
- Turn-off NIM-BIN containing Canberra modules
- Turn-on NIM-BIN containing Canberra modules
- Observe boot-up cycle on LEDs on Canberra AIM module and wait for it to complete
- Boot-up cycle is complete when “OK”-LED on Canberra AIM module comes on
- Turn-on VME crate and wait for it to reboot

2.2 MCA Spectrum Acquisition

- Control Acquisition through <Erase&Start>, <Start>, <Erase> and <Stop> Buttons in MCA window

- Spectra can be saved to a file by issuing the following command from a terminal:

Oswin% caget 27idgen1:aim_adc1 >> path/file name

File structure: generating PV • #data points • data separated by • (blanks)

Example:

```
27idgen1:aim_adc1 2048 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0
0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 1 1 0 0 0 0 0 0 0 0 1 1 2 0 0 0 0 1
0 0 0 0 0 0 0 0 0 0 0 1 0 0 0 0 0 0 0 0 0 0 1 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 1 0 1 1 0 1 0
2 2 1 1 1 1 4 2 1 3 2 4 3 3 2 4 5 1 1 4 6 5 6 2 4 7 1 7 5 4 8 10 4 6 9 6 6
8 8 11 8 12 5 13 8 13 13 22 9 9 17 12 13 12 15 10 11 18 17 15 12 20 ...
```

2.3 Energy Calibration

The assumption is that identifiable fluorescence lines of known energy are present at the detector (Am-241 radioactive source with filters, Metal (EXAFS-) foils in x-ray beam, ...).

Let a line of energy E_1 appear in channel c_1 and a line of energy E_2 appear in channel c_2 , then the energy $E(c)$ at channel c is given by

$$E(c) = (c - c_1) \frac{(E_2 - E_1)}{(c_2 - c_1)} + E_1 \quad (1)$$

The Canberra ADC is designed to convert an input pulse of 10V-amplitude to one count in the highest channel, $V(c_{\max}) = 10V$. Thus, the voltage $V(c)$ corresponding to a channel c is

$$V(c) = \frac{c}{c_{\max}} 10V \quad (2)$$


and the voltage $V(E)$ corresponding to an energy E is

$$V(E) = \frac{10V}{c_{\max}} \left((E - E_1) \frac{(c_2 - c_1)}{(E_2 - E_1)} + c_1 \right) \quad (3)$$


This voltage conversion can be used to set the single-channel analyzer (SCA)

Appendix: Sample Spectra (Am-241 source with filters)


Cu Gain = 6.00




Rb Gain = 3.50




Mo Gain = 3.50



Ag Gain = 2.10



Ba Gain = 2.20



Combined **Cu, Rb, Mo, Ag** Gain = 3.20

