Operation of the Diamond Phase Plate at Beamline 27-ID

Under construction, last updated 08/10/2016, tg

Contents

0. Principle of Operation (Bragg Transmission Phase Plate)	1
1. Hardware Setup (see Figure 1 and 2 below)	
1.1 Attocube Rotation Stage	
1.2 Kohzu Vertical Stage	
1.3 Kohzu Horizontal Stage	
1.5 IXOIIZU 11011Z011tut Diugo	/

0. Principle of Operation (Bragg Transmission Phase Plate)

By utilizing a diffraction plane rotated by 45° from the vertical around the incident synchrotron radiation beam, its natural horizontal polarization is divided into two equal in-plane (σ ') and out-of-plane (π ') components (see figure below). In a Bragg reflection these components are diffracted differently, generating a phase shift between them (Batterman, PRB **45**, 12677 (1992)). For a phase shift of π , a vertical, linear polarization results.

1. Hardware Setup (see Figure 1 and 2 below)

1.1 Attocube Rotation Stage

Back panel of ANC300 controller

- connect Attocube rotation stage to AXIS 2 output
- connect Ethernet cable on controller and on available network jack in 27-ID-B enclosure

Front panel of ANC300 controller

- connect DC-IN of ANM300 module to output of a beamline control DAC channel

Other

- connect PIN detector to input of an SR570 current amplifier > Nova1200 VFC > VME counter




Figure 1




Figure 2

Network Connection for ANC300 Controller

Network Name: anc300-27


IP #: 164.54.126.32 Mask: 255.255.255.0 Gateway: 164.54.126.1

Connect to ANC300 Controller from LINUX, using Telnet



Authorization Code: 123456

Connect to ANC300 Controller from Windows, using Putty



```
anc300-27.xray.aps.anl.gov - PuTTY
attocube ANC300 controller version 1.1.0-1304 2013-10-17 08:16
ANC150 compatibility console

Authorization code: ******
Authorization success
> stepu 2 1000
OK
> stepd 2 1000
OK
> |
```

Authorization Code: 123456

1.2 Kohzu Vertical Stage

Type: ZA10A-W2C


Motion Range: $\pm 2.5 \text{ mm}$ Lead Screw Pitch: 1 mm / rev

Motor Full Step: 0.36°/step, 1000 steps/rev (DATA1/DATA2 switch on driver: 0)
Motor Half Step: 0.18°/step, 2000 steps/rev (DATA1/DATA2 switch on driver: 1)


Vertical Travel: 0.25 mm / rev (1:4 wedge)

Sample EPICS Motor Setup

Half Step (Driver: 1)



Full Step (Driver: 0)



1.3 Kohzu Horizontal Stage


Type: XA10A-R1 Motion Range: ± 12.5 mm Lead Screw Pitch: 0.5 mm / rev

Motor Full Step: 0.36°/step, 1000 steps/rev (DATA1/DATA2 switch on driver: 0)
Motor Half Step: 0.18°/step, 2000 steps/rev (DATA1/DATA2 switch on driver: 1)

Horizontal Travel: 0.5 mm / rev

Sample EPICS Motor Setup

Half Step (Driver: 1)



Full Step (Driver: 0)

