

# Making the Most of APS Enhancements in Sector 2

Friday, 5 April 2002

Ian McNulty and Barry Lai

Experimental Facilities Division



Synchrotron Radiation Instrumentation  
Collaborative Access Team

# Micro-Techniques Group (SRI-CAT/Sector 2)

---

## MISSION


Development and application of instrumentation and techniques for high resolution x-ray microscopy.

## Personnel

9 Physicists  
2 Postdocs  
2 Engineers  
3 Technicians  
1 Grad Student

## SCOPE

Material, biological, and environmental science at dedicated beamlines in the 0.5-4 and 5-35 keV energy regions.





## Sector 2 programs

---

- **X-ray microscopy**
  - Transmission
  - Fluorescence
  - Microdiffraction**2-BM, 2-ID-B, 2-ID-D**
- **Coherent x-ray scattering****2-ID-B, 2-ID-D**
- **Deep x-ray lithography****2-BM-B**

# SRI-CAT Beamline 2-BM





# X-ray microscopy

---

- **Provides wealth of information**
  - **Transmission**                      measure electron density
  - **Fluorescence**                      measure elemental distribution
  - **Spectroscopy**                      extract chemical information
  - **Tomography**                      reveal 3D structure
  - **Diffraction**                      reveal structure, strain
- **Image the internal structure of thick or optically opaque objects**
- **No special sample preparation (staining, drying, sectioning)**
- **Non-destructive, quantitative, *in vivo*, *in situ***

# Recent applications

---

## Biological/biomedical

- Uptake of Pt anti-cancer agents in ovarian cancer cells (D. Phillips)
- Microbiology of bacteria (K. Kemner)
- Role of Fe and K in mycobacteria infection (D. Wagner)
- DNA photocleavage by  $\text{TiO}_2$  in mammalian cells (G. Woloschak)
- 3D micro-anatomy of crustacean crania (M. Fahrenstock)





## Environmental

- Elemental makeup of particles in air-borne pollutants (D. Cohen)
- Study of plant-fungi symbiotic relationship (S. Pratt)

## Materials science

- Microstructure of membrane protein crystals at room temperature (M. Caffrey)
- 3D study of electromigration voids in interconnects (Z. Levine)
- 3D fluorescence study of nuclear fuel shell (G. Ice)
- Study of antiferromagnetic domains in Cr (P. Evans)

# Fast x-ray microtomography at 2-BM






# Micro/nano focusing capabilities

---

- Zone plates provide best resolution compared to K-B optics and refractive lenses
- World-class performance of Sector 2 x-ray microscopes
- Development of key supporting instrumentation, methods
- Supporting other APS beamlines at existing and new CATs  
Instrumentation design (e.g., Nano-CAT)  
zone plates (in use by 6 CATs)
- No other beamlines (yet) optimized for these capabilities in SRI-CAT or at APS


# Highest resolution zone plates to date

---



$\delta_r = 40 \text{ nm}$ ,  $R = 49 \text{ }\mu\text{m}$ ,  $t = 0.13 \text{ }\mu\text{m}$ , Ni  
 $f = 5.7 \text{ mm}$  and  $\eta = 3.0\%$  at  $1.8 \text{ keV}$

# Improving zone plate efficiencies



**Relative performance of one zone plate vs. two stacked zone plates with outer zone width of 100 nm, at 8 keV.**

## Efficiency


**12%** @ 8 keV (0.8  $\mu\text{m}$  Au)  
**33%** @ 8 keV (1.6  $\mu\text{m}$  Ni)  
**39-45%** @ 8 keV (3-step blazed)  
**25%** @ 40 keV (2 stacked)

## Focus-to-background

**$3.6 \times 10^3$**  (3-step blazed)


**True x-ray lens!**

# 2-ID-D/E hard x-ray microscopy beamline




(only fluorescence setup shown for 2-ID-D)

## 2-ID-D hard x-ray (5-12 keV) microprobe



## 2-ID-E side-branch XFM



Downstream-outboard view (beam enters from left)

# X-ray fluorescence microscopy

---

## Intrinsic advantages compared to charged particles

- $10$ - $10^3$  times higher fluorescence cross sections
- $10$ - $10^5$  times better signal-to-background ratios
- $10^{-3}$ - $10^5$  times less radiation damage for same elemental detectability

Crucial for biological samples

- Better spatial resolution for samples with thickness  $>1\text{ }\mu\text{m}$
- Selectively excite one element versus another
- Better sensitivity for chemical state and local environment
- Simpler sample preparation

Study biological samples in their natural hydrated state





## Demonstrated capabilities

- **10-100 ppb** elemental sensitivity (**100 ppm** and **10 ppm** for EIXE and PIXE)
- **30 attogram** of trace element detection sensitivity
- Micro-XANES of Mn at **3ppm** concentration



# **X-RAY IMAGING STUDIES OF THE MYCORRHIZAL FUNGUS-PLANT SYMBIOSIS**

**90% of the world's plants, including essentially all economic crops, make use of mycorrhizal associations with fungi**





Z. Cai, K. Kemner, B. Lai, H. Lee, M. Miller,  
S. Pratt, W. Rodrigues, W. Yun




# TiO<sub>2</sub>-DNA nanocomposites in mammalian cells


- Cell is transfected with TiO<sub>2</sub>-DNA nanocomposites
- DNA targets specific chromosomal region
- TiO<sub>2</sub> photocleaves DNA strands upon illumination
- Potential use in gene therapy
- Map Ti distribution using X-ray induced K<sub>α</sub> fluorescence, to quantify success rate of TiO<sub>2</sub>-DNA transfection and visualize target
- Affinity of transfected DNA to ribosomal DNA causes nanocomposites to localize to the nucleolus



# Spatially coherent flux in the 1-4 keV region



## 2-ID-B soft x-ray (1-4 keV) scanning microscope




# Fast scan-on-the fly capability


---

- Position-based "on-the-fly" scanning now routine
- Use multiscannel scaler, advance channel with position steps
- Clock-normalization recovers exact dwell time
- Scan speed
  - 2 ms/pixel
  - 200  $\mu$ s/pixel
  - stepper stage
  - piezo stage
- Target
  - 10  $\mu$ s/pixel
  - 10<sup>9</sup> ph/s  $\Rightarrow$  1% statistics/pixel

# Scanning nanotomography of Al/W/Si chips




SIRT reconstruction of 13 STXM projections ( $\pm 69^\circ$ , 57 nm steps, 1573 eV) through two-level sample. Al interconnects (green) are joined by W vias (red). Two FIB markers are at top of image. This in-plane (side) view cannot be observed directly.




Bayesian reconstruction using same data. FIB plane to upper interconnect plane distance is  $5.4 \pm 0.5 \mu\text{m}$ .


# Coherence measurement of undulator radiation




Young's experiment (top view) at 2-ID-B beamline  
 With  $E = 1.1$  keV,  $50\text{ }\mu\text{m}$  entrance slit,  $\sim 200\text{ }\mu\text{m}$  exit slit, and Young's slit separations of 10-200  $\mu\text{m}$ .



Young's slits (1.6  $\mu\text{m}$  Au, 3  $\mu\text{m}$  wide, 10  $\mu\text{m}$  apart).



Young's fringes with 20  $\mu\text{m}$  slit separation.




Young's fringes with 50  $\mu\text{m}$  slit separation.




# Coherence function at 2-ID-B

(a) Measured horizontal degree of coherence  $|\mu_{12}|$  at 1.1 keV with monochromator entrance slit at 50  $\mu\text{m}$  and exit slit at  $\sim 200 \mu\text{m}$ . Exit slit is 8.0 m from Young's slits.  $|\mu_{12}|$  is dominated by beamline optics.



(b) Measured  $|\mu_{12}|$  with exit slit at  $\sim 100 \mu\text{m}$ .  $|\mu_{12}|$  is dominated by exit slit, producing sinc-like profile.



D. Paterson, et al., *Opt. Commun.* 195, 79 (2001)

## Future plans in Sector 2

---

- Expand biological, environmental, and materials science applications of x-ray microscopy
- Increase performance of soft and hard x-ray microscopes
  - Push resolution to 30 nm and sensitivity to 1 ppb
  - Push strain and texture analysis to 1  $\mu\text{m}$  scale
  - Improve user turnaround and ease of use
- Increase performance of fast x-ray tomography system
  - Push throughput to >1 3D data set/min
  - Push 3D resolution to 100 nm
  - Deliver turnkey operation
- Develop optics and microscopes for future x-ray sources



# We are **2D**-brilliance driven

---

## Coherent flux

Photons per time, spatial mode, and spectral bandwidth

$$F_c = B \lambda^2 / 4$$

# Increasing access to undulators in Sector 2

---

## X-ray microscopes

**Diffraction-limited optics can only use a few coherent modes**

**>100 coherent modes in horizontal, <10 coherent modes in vertical direction**

**Hard x-ray beamline (2-ID-D/E) requires white-beam aperturing in horizontal**

**Soft x-ray beamline (2-ID-B) apertures pink beam after mirror deflection**

## Maximum utilization of and access to undulator beams

**Rapid, easy ID branchline switching**

**Independent soft and hard x-ray undulators**

**Split ID beams with pick-off mirror**

**Also split ID sources?**

# Technical challenges to ID beam-splitting

---

## Pick-off mirror (M2B)

Horizontally deflecting, horizontally focusing

Max power ~ 600 W, power density ~ 0.3 W/mm<sup>2</sup>

Require <3  $\mu$ rad figure error control

High power load near downstream edge

Edge quality must not adversely affect both soft, hard x-ray beams

## Steering mirror (M3B)

Horizontally deflecting, vertically focusing

Max power ~ 20 W, power density ~ 0.3 W/mm<sup>2</sup>

Require <3  $\mu$ rad figure error control

Must be thin (~1 cm) in order to pass hard x-ray beam to 2-ID-D/E

# Current source limitations

---

- **X-ray microscopes require  $< 1$   $\mu\text{rad}$  beam stability over days**  
(few- $\mu\text{rad}$  stability currently, except after refill during top-up,  
when it may be as bad as 30  $\mu\text{rad}$  in H or V planes)
- **Most applications are brilliance-hungry, need  $10^2$  -  $10^3$  more!**  
(but, can make 10-100x better use of B with better detectors)
- **Large dead-time, dyn. range only  $\sim 10^6$  with current fill pattern**
- **Undulator access is precious!**

# Optimum solutions

---

- ID improvements offer greatest brilliance gain!
- Reduced horizontal emittance is next-best solution to more B;  
want more symmetrical source (don't reduce vert. emittance)
- Many-bunch mode offers 10-50x more dynamic range in fast  
counting experiments (scanning transmission x-ray microscopy)
- Need full-time XBPM feedback to achieve needed stability
- Higher current buys us *very little*.