

In-situ Pair Distribution Function Analysis: Structural Studies that Probe both Local and Long Range Structure

Peter J. Chupas

Materials Science Division, Argonne National Laboratory

Argonne National Laboratory



A U.S. Department of Energy
Office of Science Laboratory
Operated by The University of Chicago



Considerations

What does it take to do in-situ PDF measurements?

High Energy X-rays, Area Detectors

Why PDF?

Catalytic reactions occur on the Å scale.

Modeling Structural Data: How can we get more information from our data? As experimentalists should we consult those with expertise in theory and computation?


Applications of in-situ measurements

Phase transitions

Chemical Reactions (Catalysis)

Solid State Reactions

Area Detectors and the PDF method



Why use *Area Detectors*?


- With an area detector we can collect all data simultaneously.
- High energy X-rays are a necessity to probe large values of Q over small scattering angles (small detector areas).
- Counting statistics; averaging over a large solid angle allows adequate statistics for the contribution from Compton scattering to be accurately subtracted.

Chupas, Qiu, Hanson, Lee, Grey, Billinge *Journal of Applied Crystallography* 36 (2003) 1342.

Nickel Collected with an Image Plate



- 1 second counting
- A conventional measurement takes ~8 hours... we need time resolution better than ~1 minute for in-situ studies

QuickTime™ and a
TIFF (LZW) decompressor
are needed to see this picture.



In-situ Powder Diffraction



- AlF_3 undergoes a phase transition from a rhombohedrally distorted ReO_3 structure to a cubic structure at $\sim 470^\circ\text{C}$.



OpenStax and a
CC BY-SA license
are needed to use this image.





The Structure of $\alpha\text{-AlF}_3$





Rotations About the 3-fold Axis

The rhombohedral distortion in $\alpha\text{-AlF}_3$ involves rotations of rigid octahedra



Dynamics in AlF_3 at High Temperature



PDF is a measure of the “instantaneous” structure, whereas Rietveld yields the time averaged position of atoms




In-situ Pair Distribution Function Analysis

QuickTime™ and a
PDF (Acrobat) document
are needed to see this picture.

Chupas, Chaudhuri, Hanson et al. *Journal of the American Chemical Society* 126 (2004) 1342.



Using a Combination of Analysis Methods




$$R_{\text{rietveld}} = 5.55 \%$$

$$R_{\text{PDF}} = 21.51 \%$$



$$R_{\text{rietveld}} = 5.80 \%$$

$$R_{\text{PDF}} = 18.60 \%$$



Evidence for Motion of Rigid AlF_6 Octahedra


The experimental results are indicative of a dynamic model where the octahedra are rigid.




QuickTime™ and a
TIFF (LZW) decompressor
are needed to see this picture.

Snapshots from the Molecular Dynamics Simulations

- Molecular Dynamics simulations (Santanu Chaudhuri, Mark Wilson, Paul Madden)
- MD is ideally suited for comparison with total scattering methods.



500 K (Rhom)



900 K Cubic

Chaudhuri, Chupas, Wilson, Madden, *Grey Journal of Physical Chemistry* 108 (2004) 3437.



Snapshots from the MD

- Molecular Dynamics simulations (Santanu Chaudhuri, Mark Wilson, Paul Madden)

Chaudhuri, Chupas, Wilson, Madden, *Grey Journal of Physical Chemistry* 108 (2004) 3437.



Partial Radial Distribution Functions from MD

Al-F

MD allows the determination of partial pair correlations

- The Al-F bond distance does not show significant changes as a function of temperature
- However, broadening of the Al-F correlation is clearly observed.

QuickTime™ and a
TIFF (LZW) decompressor
are needed to see this picture.

Al-Al

- The Al-Al gradually increase as a function of temperature

Correlation Functions from MD


$$\mathbf{r}_c = (\mathbf{r}_k - \mathbf{r}_l)/2 + \mathbf{r}_l$$

$$\delta = \mathbf{r}_i - \mathbf{r}_c$$


$$C_{\text{disp}}(t) = \langle |\delta_i(t)| |\delta_i(0)| \rangle$$

$$C_{\text{angle}}(t) = \langle |\hat{\delta}_i(t) \cdot \hat{\delta}_i(0)| \rangle$$



QuickTime™ and a TIFF (Uncompressed) decompressor are needed to see this picture.





AlF₃: An example of an Industrial Important Heterogeneous Catalyst



nano-AlF₃ from PDF and MD




Catalytically inactive $\alpha\text{-AlF}_3$



Catalytically inactive $\beta\text{-AlF}_3$

Studying Local Structure During Chemical Reactions


Utilize a reaction cell designed to accurately control temperature and control the flow of gases over a sample.



Ceria as a Catalyst

Ceria is a major component of “Three Way” Catalysts: Which Simultaneously treat the reducing pollutants CO and C_xH_y , and the oxidizing pollutant NO_x

Oxidation: $2 CO + O_2 \rightarrow CO_2$





Reduction $2 CO + 2 NO \rightarrow 2CO_2 + N_2$



Ceria acts as an oxygen reservoir to stabilize the air/fuel ratio

Fluorite structure with variable composition





In-situ Reduction of Ceria


No defect ordered phases identified;

$\epsilon\text{-CeO}_{1.770}$




$\delta\text{-CeO}_{1.817}$



Studying Local Structure During Chemical Reactions





Frenkel Defects in nano-CeO₂



Mamontov and Egami *Journal of Physics and Chemistry of Solids* 61 (2000) 1345.

Nano-CeO₂: Changing Oxygen Content

Oxygen is removed from the 01 site in nano-Ceria on Reduction



Faster Area Detectors



For PDF measurements we need to be able to average over a large solid angle to obtain adequate counting statistics.

Faster Readout: 7-30 Hz

Image Plate ~1min.



CCD ~ 4 sec.

Studying Local Structure During Chemical Reactions


- Potential Hydrogen Storage Material



Diffraction experiments with 1 second resolution necessary to isolate amorphous phase



Studying Local Structure During Chemical Reactions



Final Thoughts

- The simple conclusion: In-situ experiments are critical to correlate structure/property relationships.

What is needed to improve and expand the application of in-situ studies?

Instrumentation

High Energy X-rays, Detectors

Advances in Modeling

Incorporating advances in modeling techniques

Ancillary Equipment

Optimized control of the sample environment

Acknowledgements

•Clare Grey	SUNY-Stony Brook
•Santanu Chaudhuri	SUNY-Stony Brook (currently BNL)
•Paul Madden	Oxford University
•Mark Wilson	University College London
•Jon Hanson	Brookhaven National Laboratory
•Peter Lee	APS, Argonne National Laboratory
•Sarvjit Shastri	APS, Argonne National Laboratory
•Yan Gao	General Electric
•Kenneth Kump	General Electric
•German Vera	General Electric
•C. Shawn Rogers	General Electric
•Simon Billinge	Michigan State University
•Xiangyun Qiu	Michigan State University

•Work supported by Department of Energy, Office of Basic Energy Sciences. Work at the APS was supported by U.S. Department of Energy, Office of Science under contract W-31-109-ENG-38

