

# ***The Present and Future Use of High-Energy X-rays for Industrial Materials Research***

*Yan Gao  
GE Global Research Center  
Niskayuna, NY 12309*

*Workshop on "Science with High-Energy X-rays"  
August 9, 2004, Advanced Photon Source, ANL*

# Acknowledgement



*Beamlines XOR 1-ID-C, 5-BM-D*

*Ulrich Lienert (APS)*

*Jon Almer (APS)*

*Peter Lee (APS)*

*Dean Haeffner (APS)*

*Peter Chupas (ANL)*

*Qing Ma (DND-CAT)*

*Beamlines X17B1, X15A*

*Zhong Zhong (NSLS)*

*Bill Carter (GE)*

*Jim Ruud (GE)*

*Tom Angeliu (GE)*

*Kan Kump (GEMS)*



# Outline

## ➤ *GE and GE Global Research*

### ➤ *The Present*

- *Residual stress measurement*
- *Characterization of TBC*
- *High throughput XRD and SAXS*
- *Hg XRF at 83 keV*
- *Pr EXAFS at 42 keV*
- *Other applications with HE X-rays*

### ➤ *The Future*

- *Cutting-edge capability*
- *Advanced characterization*
- *A friendly user facility*

# ***GE and GE Global Research***



***Niskayuna, NY***



***India***



***China***



***Germany***

imagination at work



# ***GE Global Research: Hub for innovation***

***Then...***




***Discovery of synchrotron radiation at GE Research Center (1947)***

***And Now...***

***Cutting-edge research***

- ***Nanotechnology & functional materials***
- ***Hydrogen storage materials***
- ***Solid Oxide Fuel Cell***
- ***Photovoltaics***




imagination at work






# Why is GE interested in HE X-rays

- *Unique capability (penetrating power. intensity)*
- *Superior data quality (S/N, angular resolution)*
- *Productivity (simpler sample prep, fast data collection)*



*Research on many metals, alloys, and ceramics*

# Residual stress measurement



Obtain depth profile by layer removal




## In-house X-ray source

- Large beam footprint
- Low intensity for high-angle peaks
- Low accuracy
- Layer-removal for depth profile

## Synchrotron X-rays

- Small beam size
- High intensity
- High accuracy
- Non-destructive with HE X-rays




Obtain depth profile by moving sample

*Residual Stress Determination Is Very Important for Industry*


# Non-destructive residual stress measurement

- Shot-peening effect was measured as a function of depth
- Sample stage was used to bring more grains to diffraction

## Obtaining triaxial strain tensor



Diffraction data were collected at various  $\chi$ ,  $\phi$  and  $t$  (depth up to 1.3 mm)  
352 images were taken in 6 hours per sample in automated operation




XOR 1-ID-C (U. Lienert)



# Non-destructive residual stress measurement


- Triaxial stress components determined
- Affected depth up to 400  $\mu\text{m}$ .



# High-energy X-ray powder diffraction

*Combined with MarCCD and an auto sampler, high throughput measurements can be performed.*

*Sagittal focusing monochromator*



- Providing  $10^{11}$  ph/s at 67 keV by focusing horizontal beam from 40 mm to 0.5 mm, which is an flux increase by 2 order of magnitude.
- Vertical divergence is between 10-30 micro-radians good angular resolution

*X17B1 (Z. Zhong)*

# Applications to Thermal Barrier Coatings (TBC)




➔ YS  
Z

➔ Bond coat


➔ Superalloy substrate

- **Polymorphs: tetragonal, cubic and monoclinic**
- **Separation of tetragonal and cubic peaks**
- **Determination of lattice parameters and  $c/a'$ : transformable (t) and non-transformable (t')**
- **Texture: difficult with conventional XRD.**



# Applications to Thermal Barrier Coatings (TBC)

*TBC analyses may involve large number of measurements*





or

## Performance DOE

Time	Temperature		
	x	x	x
	x	x	x
	x	x	x
	x	x	x
	x	x	x


- Tetra = 49.9 wt%
- Cubic = 47.9 wt%
- Mono = 2.2 wt%
- Tetra  $c/a' = 1.0154 \rightarrow t'$
- Tetra  $c/a' \rightarrow 4.4 \text{ mol\% } YO_{1.5}$
- Tetra  $c/a' \rightarrow \text{thermal history}$
- Cubic  $a \rightarrow 14.0 \text{ mol\% } YO_{1.5}$
- Peak width  $\rightarrow \text{micro-strain}$
- Peak position  $\rightarrow \text{macro-strain}$




# Applications to cast ODS alloys




## Oxide Dispersion Strengthened Alloys



## Conventional analyses




- XRD phase analysis
- TEM sample: from a tiny area
- SAXS sample: thin and small



TEM micrograph




# Applications to cast ODS alloys




## Transmission HE-XRD

- Phase identification
- Oxide dispersion in macro scale



## HE-SAXS

- Oxide dispersion in micro scale
- Oxide size and size distribution




# Applications to cast ODS alloys




## HE-SAXS

- Oxide dispersion in micro scale
- Oxide size and size distribution



*Igor + SAS macro (J. Ilavsky)*



*Size and size distribution*

*XOR 1-ID-C (J. Almer)*

## High-energy X-ray fluorescence

### Example: Hg K-edge

**IA** **IIA** **IIIA** **IVA** **VA** **VIA** **VIIA** **0**

**1 H** 1.008 **2 He** 4.003

**3 Li** 6.941 **4 Be** 9.012 **5 B** 10.81 **6 C** 12.01 **7 N** 14.01 **8 O** 16.00 **9 F** 19.00 **10 Ne** 20.18

**11 Na** 22.99 **12 Mg** 24.31 **13 Al** 26.98 **14 Si** 28.09 **15 P** 30.97 **16 S** 32.06 **17 Cl** 35.45 **18 Ar** 39.95

**19 K** 39.10 **20 Ca** 40.08 **21 Sc** 44.96 **22 Ti** 47.90 **23 V** 50.94 **24 Cr** 52.00 **25 Mn** 54.94 **26 Fe** 55.85 **27 Co** 58.93 **28 Ni** 58.70 **29 Cu** 63.55 **30 Zn** 65.38 **31 Ga** 69.72 **32 Ge** 72.59 **33 As** 74.92 **34 Se** 78.96 **35 Br** 79.90 **36 Kr** 83.80

**37 Rb** 85.47 **38 Sr** 87.62 **39 Y** 88.91 **40 Zr** 91.22 **41 Nb** 92.91 **42 Mo** 95.94 **43 Tc** (98) **44 Ru** 101.1 **45 Rh** 102.9 **46 Pd** 106.4 **47 Ag** 107.9 **48 Cd** 112.4 **49 In** 114.8 **50 Sn** 118.7 **51 Sb** 121.8 **52 Te** 127.6 **53 I** 126.9 **54 Xe** 131.3

**55 Cs** 132.9 **56 Ba** 137.3 **57 La** 138.9 **72 Hf** 178.5 **73 Ta** 180.9 **74 W** 183.9 **75 Re** 186.2 **76 Os** 190.2 **77 Ir** 192.2 **78 Pt** 195.1 **79 Au** 197.0 **80 Hg** 200.6 **81 Tl** 204.4 **82 Pb** 207.2 **83 Bi** 209.0 **84 Po** (209) **85 At** (210) **86 Rn** (222)

**87 Fr** (223) **88 Ra** (226.0) **89 Ac** (227) **104 Rf** **105 Ha** **106 Unh** **107 Uns** **108** **109 Une**

**90 Th** 232.0 **91 Pa** (231) **92 U** 238.0 **93 Np** (244) **94 Pu** (242) **95 Am** (243) **96 Cm** (247) **97 Bk** (247) **98 Cf** (251) **99 Es** (252) **100 Fm** (257) **101 Md** (258) **102 No** (259) **103 Lr** (260)

**36.0** **90.5**




- Many commonly used elements, including RE elements
  - Excitation beyond the energy of in-house XRF
  - Greater fluorescence yield at K-edge
  - Useful for non-destructive detection
- Non-destructive
  - Hg vapor pressure vs. temperature





# Non-destructive detection of Hg vapor in F-lamps

Setup at NSLS X26A and X17B1




The lamp was wrapped with heating tape for temperature-dependent measurement



# X-ray absorption spectroscopy using HE X-rays

**Motivation: Understand the role of Pr doping in Quantum Splitting Ph**




- Pr L-edge not possible due to La
- Pr K-edge not possible at

	38.91	91.22	92.91	
57	La	72	Hf	73
138.9		178.5		180.9
89	Ac	104	Rf	105
(227)				106
				Unh
				U


58	59	60	
Ce	Pr	Nd	F
	140.9	144.2	
	91	92	



**K-edge data from APS 5-BM-D**





# Quantitative phase analysis using HE X-rays



Objective: Quantify W, WC and W<sub>2</sub>C

*Tungsten absorption is too severe for in-house conventional X-rays*

Solution: High-energy XRD at 67 keV!



# Non-destructive XRD using HE X-rays



## Patents

United States Patent 6,132,823  
Giu October 17, 2000

## Superconducting heat transfer medium


### Abstract

A superconducting heat transfer medium that has three basic layers, the first layer being various combinations of sodium, beryllium, a metal such as manganese or aluminum, calcium, boron and dichromate radical; the second layer formed over the first layer and being various combinations of cobalt, manganese, beryllium, strontium, rhodium, copper, beta-titanium calcium, a metal such as manganese or aluminum and the dichromate radical.




*How does it work and why?*

*The pipe is claimed to have several layers, and work only when it is sealed, therefore **HE-XRD** is the chosen technique to investigate the interior chemistry and crystal structure.*



*67 KeV X-ray in transmission mode*



# *The Future*

## *Cutting-edge capability*

- *A turnkey facility for MicroXRD and microXRF*
- *Fast time-resolved in-situ diffraction*

## *Advanced characterization*

- *Non-destructive residual stress*
- *High-throughput materials screening (XRD, XRF and SAXS)*

## *A friendly user (including industrial users) facility*


- *Dedicated instrumentation for frequently used techniques*
- *Quick access and/or remote access*
- *Commercialized analytical services*



# Microdiffraction with high spatial resolution


**Cross-section of a SOFC part:**

**Consisting cathode, anode, electrolytes, interconnect, seal glass,**




*Optical micrograph*

**While elemental information may be obtained by SEM-EDS, it's very important to obtain crystal structure information from region of interest as marked.**




# A turnkey microdiffraction station: aim and shoot




Monochromatic beam for polycrystalline area



White beam for single crystal grains



A wavelength dispersive XRF detector can be added



Crystal structure




# Fast time-resolved XRD



*Time resolution can be essential for mechanistic understanding!*



*~70 sec, 21 measurements*



*XOR 1-ID-C using GE 2D detector (P. Chupas and P. Lee)*




# Fast time-resolved XRD with GE Detector

*GE detector at work at XOR 1-ID-C*



- Area: 41 cm x 41 cm
- Pixel size: 200  $\mu\text{m}$
- Readout: 41 ms for 2k x 2k (Angio)
- Dynamic range: 14 bits




Contact: Ken Kump, GE Medical Systems, Ph: 262-548-4549

imagination at work



# Advanced Characterization

## High throughput XRD measurements




*With automated sample stage, 2D detector, high-energy X-rays, superb synchrotron intensity, and dedicated data analysis software, large number of samples can be preloaded and measured unattended or remotely. High-energy X-rays in transmission mode is particularly useful for many inorganic or metallic materials.*


*Same approach, combined with micro-focused beam, can also be used for automated diffraction mapping with monochromatic beam, or elemental mapping with white beam.*

# Advanced Characterization


- **Residual stress** and **plastic deformation** are very important for industrial applications
- Actual samples involves complex geometry, and small beam and high intensity are essential for obtaining accurate data
- **Non-destructive with HE X-rays**
- Consider dedicated instrumentation and commercialize the analyses through a third party.



## Residual stress



## Plastic deformation



# The Future

*A friendly user (including industrial users) facility*

- *Dedicated instrumentation for frequently used techniques*

*Such as powder diffraction with 2D detector for normal or high throughput applications*

- *Quick access and/or remote access*

*Linked with dedicated instrument to minimize setup time; web-based remote access for users running experiment from home institution*

- *Commercialized analytical services*

*Powder diffraction and residual stress measurement may be two key areas to promote fee-based analytical services*