

# Site-selective magnetization measurements using nuclear resonant scattering

**C. L'abbé**

Instituut voor Kern- en Stralingsfysica, Katholieke Universiteit Leuven, Belgium  
[caroline.labbe@fys.kuleuven.ac.be](mailto:caroline.labbe@fys.kuleuven.ac.be)


Outline :

- I. Introduction
- II. Nuclear resonant scattering
- III. Circularly polarized radiation
- IV. Interlayer coupling in Fe/Cr multilayers
- V. Conclusions



Site-selective magnetization measurements :

- XMCD

- element-specific scattering
- study different materials independently



- Nuclear resonant scattering




- isotope-selective scattering
- study specific sites within the material separately

## II. Nuclear resonant scattering

Nuclear resonant scattering of synchrotron radiation :

- **isotope-selective** technique : only scattering on the resonant nuclei ( $^{57}\text{Fe}$ )
- measures the local hyperfine interaction in the sample
  - isomer shift                      ~ chemical environment of probe nuclei
  - electric field gradient              ~ lattice symmetry around the probe nuclei
  - **magnetic hyperfine field**              ~ **magnetization properties**
- measures the resonantly scattered synchrotron radiation **versus time**




# Information in time spectra


The time spectrum is

- sensitive to the magnitude of the magnetization vector

$$M = 1.76 \text{ J/Tm}^3$$



$$M = 0.88 \text{ J/Tm}^3$$





⇒ beat frequency  $\sim$  magnitude of  $M$

# Information in time spectra

The time spectrum is





- sensitive to the magnitude of the magnetization vector
- sensitive to the direction of the magnetization vector



# Information in time spectra

The time spectrum is

- sensitive to the magnitude of the magnetization vector
- sensitive to the direction of the magnetization vector
- restriction :



⇒ two opposite directions of  $M$  give exactly the same time spectrum !

The time spectrum is





- sensitive to the magnitude of the magnetization vector
- sensitive to the direction of the magnetization vector
- restriction : **NOT sensitive** to the sign of the magnetization vector

Explanation :

because the incident radiation is **linearly polarized**  
the scattering process is **not sensitive to the sign of M**

# III. Circularly polarized radiation

In order to measure also the **sign of M**  
one has to use **circularly polarized** radiation




⇒ two opposite directions of M give clearly distinct time spectra !




# IV. Interlayer coupling in Fe/Cr multilayers

Fe/Cr multilayers :




Depending on the Cr layer thickness, the Fe magnetization vectors will align :

- under  $0^\circ$  or  $180^\circ$  : **bilinear coupling**




- under  $90^\circ$  : **biquadratic coupling**




In order to study the interlayer coupling in detail :

- simplify the multilayer system :



- quintalayer grown on MgO(100)

standard magnetization measurement :




⇒ measures the response of all Fe layers together

How can one measure a magnetization curve of 1 layer selectively ?

In order to study the interlayer coupling in detail :

- simplify the multilayer system :
- measure magnetization of each Fe layer independently




Use the isotope-selectivity of nuclear resonant scattering :




- buried  $^{57}\text{Fe}$  layer
- quintalayer grown on thick  $\text{MgO}(100)$

Measurement yields the magnetization vector of the central Fe layer


# Experiment at APS beamline 3-ID



- Sample grown with MBE on MgO(100) :




## Nuclear resonant magnetometry




- ⇒ isotope-selective magnetization curve
- ⇒ reflects only the central Fe layer

- Sample grown with sputtering on MgO(100) :



⇒ at zero field, the central magnetization vector is NOT antiparallel !!


Retrieve quantitative values for **coupling angle** :



$\phi$  : angle of outer magnetization vectors

$\theta$  : angle of central magnetization vector


nuclear resonant magnetometry :



central Fe layer

$$M/M_S = \cos \theta$$


standard magnetometry :



all Fe layers

$$M/M_S = (2\cos \phi + \cos \theta)/3$$


Retrieve quantitative values for **coupling angle** :



$\phi$  : angle of outer magnetization vectors

$\theta$  : angle of central magnetization vector

Combining the **nuclear resonant magnetization data** and the **standard magnetization data** detailed information on the interlayer coupling can be obtained



at zero field :  $|\theta - \phi| = 162^\circ \pm 4^\circ$

⇒ non-collinear coupling !!




## V. Conclusions




- Nuclear resonant scattering **with circularly polarized radiation** permits to retrieve **detailed magnetic information**
- We measured a layer-selective magnetization curve in  $[\text{Fe}(5.0\text{nm})/\text{Cr}(1.1\text{nm})]_3$  and found
  - bilinear coupling for MBE-grown samples
  - **non-collinear** coupling for sputtered samples

# Perspectives : other applications

- **depth-selective** measurement of **spinrotation** in exchange-coupled bilayers :






⇒ alignment of ferromagnet by antiferromagnet



# Perspectives : other applications

- depth-selective measurement of spinrotation in exchange-coupled bilayers :
- measurement of **nanoscale islands** :



⇒ perpendicular spin orientation in Fe islands below 100 K

## Collaborators :

Instituut voor Kern- en Stralingsfysica, University of Leuven, Belgium

J. Meersschant

High-Resolution X-ray Scattering, Advanced Photon Source, Argonne National Lab.

W. Sturhahn

T.S. Toellner

E.E. Alp

Magnetic Thin Films, Materials Science Division, Argonne National Laboratory

J.S. Jiang

S.D. Bader

Work at Leuven was supported by the Fund for Scientific Research Flanders and the Inter-University Attraction Pole IUAP P5/1

Work at Argonne and the use of the APS was supported by U.S. DOE, BES Office of Science, under Contract No. W-31-109-ENG-38