

# Surface analysis of biological materials using ToF-SIMS and laser post-ionisation

Nick Lockyer & John Vickerman

- Introduction
  - ToF-SIMS
  - Laser Post-Ionisation
- Analysis of biological materials
  - Challenges
  - Human Prostate Cancer Cells
  - Laser PI & the FEL
- Summary & Conclusions


Surface  
Analysis  
Research  
Centre



# Secondary Ion Mass Spectrometry (SIMS)


## *Laser Post-ionisation (Sputtered Neutral Mass Spectrometry)*

- *Primary* ion beam sputters atomic and molecular species from *top monolayer* of a solid surface
- 'Static SIMS' (<1% impacted) – detailed molecular characterisation
- *Secondary* ions are extracted and analysed by mass spectrometry





# Retrospective Imaging - Chemical Microscopy

- Digitally scan focused ion beam over the surface
- FULL mass spectrum at each pixel.**




# ToF-SIMS - Instrumentation



# Laser Post-Ionisation - Mechanisms

- Single photon ionisation (SPI)
- Multiphoton ionisation (MPI)
  - resonant (REMPI)
  - non-resonant (NRMPI)
- Field Ionisation (FI)
  - high-intensity MPI
  - tunnel ionisation (TI)
  - barrier-suppression (BSI)


# Laser Post-Ionisation



# Bioanalysis – Challenges

- Sample preparation
  - Vacuum compatibility of biological cells
- Data interpretation
  - Complex mass spectra
  - Fragmentation
- Sensitivity
  - Low concentrations
  - Low secondary ion yields
  - Fragmentation

# Sample Preparation – Fast Freezing & Freeze Fracture




- *in-situ* fracture (frozen hydrated surface for analysis)
- or *ex-situ* fracture & freeze dry




# ToF-SIMS Characterisation of Human Prostate Cancer cells

- Gleason *histopathological* Grading System
  - based on glandular architecture


Grade 2




Grade 3



Grade 4



Grade 5




*Disease progression*

- Need a **chemically based** diagnostic tool to study basis of carcinogenesis


# Data interpretation – Univariate Approach

- ToF-SIMS images of freeze-dried PC-3 cancer cells

**Intact** 150  $\mu\text{m}$  f.o.v




**Fractured** 100  $\mu\text{m}$  f.o.v




# Data interpretation – Multivariate Approach

- Principal Component Analysis of cell lines

ToF-SIMS differentiates non malignant and malignant cancer cells derived from different metastatic sites




# Loadings plot for principal components 1-3 from CaP cell lines



- Helps interpret the scores plot
- Provides clues to the biochemistry underlying metastatic potential



# Sensitivity – Molecular imaging



Pixel size	Pixel area	Molecules per pixel *	Atoms per pixel
10 $\mu\text{m}$ x 10 $\mu\text{m}$	$10^{-6} \text{ cm}^2$	$4 \times 10^8$	$2.5 \times 10^9$
1 $\mu\text{m}$ x 1 $\mu\text{m}$	$10^{-8} \text{ cm}^2$	$4 \times 10^6$	$2.5 \times 10^7$
500 nm x 500 nm	$2.5 \times 10^{-9} \text{ cm}^2$	$1 \times 10^6$	$6.25 \times 10^6$
100 nm x 100 nm	$1 \times 10^{-10} \text{ cm}^2$	40 000	$2.5 \times 10^5$
200 Å x 200 Å	$4 \times 10^{-12} \text{ cm}^2$	1600	10 000

\* assuming a molecular area of 5 Å x 5 Å

- ToF-SIMS**

- Static conditions (< 1% removed)  $\Rightarrow 10^4$  molecules for analysis
- Secondary ion yield  $\sim 10^{-5}$ - $10^{-3} \Rightarrow \sim 0.1$ - $10$  ions/pixel

# Sensitivity – increased ion yield


$$I_s^{+/-} = I_p Y_{chem} \alpha^{+/-} \theta_{chem} T$$

- Secondary ion yield of biomolecules (1k-10kDa) is very low  $< 10^{-5}$
- 1990s: Polyatomics increase ions yields from sputtered organics
  - Texas, large hydrocarbons,  $C_{60}^+$
  - Orsay,  $Au_n^+$
  - Munster, Idaho,  $SF_5^+$
- 2001: UMIST/Ionoptika Ltd develop practical ToF-SIMS ion beam sources based on  $Au_n^+$  and  $C_{60}^+$


- Non-linear yield increases for polyatomics

PET, 10 keV, + SSIMS

■ Au<sup>+</sup> ■ Au<sup>3+</sup> ■ C<sub>60</sub><sup>+</sup>




- Non-linear yield increases for polyatomics





- High-mass & polyatomic primary beams increase secondary ion yield



D. Weibel *et al.* *Anal. Chem.* 75 (2003) 1754

# Mechanism of Secondary Ion Yield Enhancement with polyatomics?

- Impact energy spread over wider area.
- Impact energy dissipated very close to surface - greatest effect with  $C_{60}$ .
  - => Multiple sputtering events
  - => Non-linear enhancement
- What is the mechanism of emission?
- What is the effect on the internal energy of sputtered molecules?
- What about the neutral yield?

Threshold ionisation with tuneable VUV?


# Sensitivity - Post-ionisation

$$I_s^0 = I_p Y_{chem} \alpha^0 \theta_{chem} T$$

- Neutral yields typically orders of magnitude higher than ion yields
- Post-ionisation efficiency determined by:
  - photoionisation efficiency  
 $\leq 100\%$  (atoms)  $\leq 1-10\%$  (organics) (SIMS  $\alpha^{+/-} \leq 10^{-3}$ )
    - Laser parameters: intensity, wavelength, pulse duration
    - Internal energy of desorbed species
  - temporal/spatial overlap of post-ionisation volume with sputter plume



# Harmonic Generation of VUV photons for SPI



$$I(3\omega) \propto N^2 |\chi^{(3)}(\omega)|^2 I(\omega)^3 F(L, b, \Delta k)$$

- Conversion efficiency ~0.01%, **10<sup>12</sup> photons/pulse** at 118 nm

A.H. Kung *et al.* *Appl. Phys. Lett.* **22** (1973) 301 & **28** (1976) 239

U. Schühle *et al.* *J. Am. Chem. Soc.* **110** (1988) 2323 (1988)

# Molecular post-ionisation - Fragmentation

## Multiphoton




- Selectivity based on 'resonance-enhancement'
  - Sputtered molecules occupy many vibrational modes
- Tuneable fragmentation
  - determined by laser parameters:  $I$ ,  $\lambda$ ,  $\tau$

## Single photon

- Non selective ionisation
- 'Soft' ionisation
- Tuneable VUV (e.g. from the FEL)
  - Threshold ionisation (increased selectivity)
  - Tuneable fragmentation




# Post-ionisation - Fragmentation

- Internal energy effects



# Post-ionisation - Fragmentation

- Ultrafast (ps-fs) pulses reduce MPI fragmentation





# Summary - Biomolecular ToF-SIMS

- Cellular samples can be successfully prepared for ToF-SIMS/PI analysis using cryo-techniques
- Chemometric data analysis aids in the interpretation of complex data sets from biological samples
- Sensitivity is still a key issue in obtaining useful molecular information at the subcellular level
  - high-mass/polyatomics primary beams help
  - **optimised post-ionisation would also help!!!**

# Future Directions

- Studies of fundamental cellular processes
- Providing information on the biochemical basis for disease progression at the cellular level
- Rapid identification of pathogenic microbes directly on food-stuffs, forensic samples etc.

# Summary – Biomolecular Laser Post-ionisation

- ✓ Decoupled desorption & ionisation
  - Independent optimisation of each step
- ✓ MPI - Some selectivity
  - multiple chromophores, uncontrolled photon absorption  
⇒ fragmentation
  - Ultrafast laser pulses reduce photofragmentation
- ✓ SPI - Non selective
  - <1% ionisation efficiency,  $10^{12}$  photon/pulse,  $10^5$  W/cm<sup>2</sup>
- ✓ Less complex mass spectra (?)
- ✓ Greater sensitivity (?)
  - Sub femtomol detection limit
- ✗ Extra cost & complexity

# Applications of Free Electron Laser in Biological SIMS/PI

- Spectroscopic investigation of sputtered ions and neutrals
  - Insight into mechanisms of sputtering regimes
- More control over ionisation and fragmentation of sputtered biomolecules
  - MPI with tuneable UV femtosecond pulses
  - SPI with tuneable VUV
- High sensitivity post-ionisation
  - $>10^{12}$  photons/pulse
  - 2-photon absorption may be a problem for SPI  $\alpha^0 > 0.3$

# Conclusion

Laser post-ionisation is a valuable complementary technique to ToF-SIMS offering many advantages in both applications and basic research.

# Acknowledgements

## ORGANIZING COMMITTEE FOR INVITATION

**FUNDING:** UK Engineering and Physical Sciences Research Council (EPSRC)

**INSTRUMENT DEVELOPMENT:** Prof. Nick Winograd, Penn State Univ.  
Kore Technology Ltd  
Ionoptika Ltd

**POLYATOMICS:** Dr. Daniel Weibel Nathan Davies

**CELL WORK:** Ehsan Gazi Emrys Jones  
Dr. Harald Jungnickel Dr. Ben Cliff  
Dr. Peter Gardner Prof. John Dwyer



## Christies Hospital / Paterson Institute

Mr. Noel Clarke  
Dr. Jonathan Shanks

Dr. Mick Brown  
Clare Hart