
The Deep Ultra-Violet Free Electron Laser

**X.J. Wang for DUV-FEL Team
National Synchrotron Light Source**


*Presented at the ALFF Workshop
October 30, 2003*

Outline


- Introduction
- DUV-FEL Facility
- Recent Experimental Results - SASE and HGHG
- DUV- FEL user science program - ion pair imaging experiment
- Future upgrade and Summary

Free Electron Laser Configurations

OSCILLATOR




SINGLE PASS FEL



Challenges:

1. Oscillator produces fully coherent output; but there is no High-quality mirrors in UV and X-ray range.
2. SASE output covers full spectrum, but it is only transverse coherent.
3. HGHG is capable of producing fully coherent output, seed laser limits its spectra coverage.

High Gain Harmonic Generation (HGHG) Principle



HGHG has the following advantages:

- Longitudinally fully coherent
- Narrower bandwidth than SASE
- Larger ratio of output/spontaneous radiation
- Central wavelength is stable
- Pulse length is short & controllable (20 fs)
- Output fluctuations can be reduced




DUV - FEL Facility

- A dedicated platform for the development of single-pass FEL Science and Technology, and its applications.
- The only short wavelength FEL project based on Laser Seeded High Gain Harmonic Generation
- Many collaborative relationships: BNL-ATF, APS, SLAC-LCLS, TJNAF-FEL, Duke, U. Maryland, TESLA


DUV - FEL Facility



- High-Brightness femto-second electron beam.
- Synchronized femto-second laser beam.
- Sophisticated electron beam and laser instrumentation
- HGHG FEL:

Electron Energy (MeV)	176		Peak Current (A)	300
NISUS Period λ_U (cm)	3.89		NISUS length (m)	10
Seed laser λ_S (nm)	800		Seed laser pulse length (FWHM) (ps)	0.1 - 6
Energy/pulse at 266 nm (μJ)	100		Energy/pulse at the 3rd harmonic 89 nm (μJ)	~ 1.0
HGHG pulse length (FWHM) (ps)	1 – 0.5		HGHG spectrum width (%)	0.1
Spot size at 266 nm (rms, μm)	250		Spot size at 89 nm (rms, μm)	150



Brookhaven Science Associates
U.S. Department of Energy






HGHG Diagnostics
Station

Undulator and Electron Beam Parameters


Period	3.89 cm
Number of periods (length)	256 (10 m)
Peak field	0.31 T
Betatron wavelength (at 140 MeV)	20 m
Electron beam size, RMS (4 mm mrad)	250 μm

Energy	Up to 200 MeV
Charge	300 pC
Normalized emittance	4 mm*mrad
Compressed bunch length, RMS	0.3-0.6 ps
Energy spread, RMS	0.3 %

SASE Signal at DUV-FEL February 2002




266 nm SASE Signal along NISUS Wiggler





$L_G = 0.66$ m
 $I_o = 550$ Amp
 $\epsilon_n = 3$ mm-mrad
 $Q = 300$ pC

266 nm HGHG Power vs. Distance




Average output: 100 μJ , 10% fluctuation

Spectrum of HGHG and SASE at 266 nm




HGFG Pulse Length Measurements




HGHG Pulse Length Measurements

Two-photon absorption pump probe
autocorrelation traces



- Pulse length is 0.63 ps
- Seed laser 1.8 MW
- Saturation at the end of wiggler







- Pulse length is 0.81 ps
- Seed laser 80 MW
- Saturation after 5 m

HGHH Pulse Length Control by Seed Laser



HGHH pulse was measured using two-photon absorption auto-correlator for seed laser 6 ps to 150 fs. Slippage could lead to 100 to 200 fs HGHH pulse lengthening. Possible measurement errors are:

- Resolution of micrometer (70 fs).
- Jitters .
- GVD Dispersion (70 fs/nm).





Intensity Fluctuation of SASE and HGHG

Unsaturated
SASE



HGHG




Harmonic Generation


$$P_2 = P_3 \left(\frac{K}{\gamma k_u \sigma_x} \right)^2 \left(\frac{K_2}{K_3} \right)^2 \left(\frac{b_2}{b_3} \right)^2$$

$$\left(\frac{P_3^{NL}}{\rho P_{\text{beam}}} \right) \approx |H_0|^2 \frac{16 w_{1r}^3}{w_{3r}} \left(\frac{P_1}{\rho P_{\text{beam}}} \right)^3 \sim 1 \%$$


undulator
entrance







half-way
saturation



full
saturation







Scaled to 100 μJ for 266 nm, assumed mono has equal efficiency at 89 and 266 nm, and a factor of 2 greater at 133 nm.

Other Capabilities at the DUV-FEL Facility

- Coherent THz Radiation
- Femto-second pulse radiolysis

Possible Future development:

- Ultra-fast plasma X-ray source.
- Femto-second electron diffraction.



XUV-FEL Scientific Opportunities

Near-term

Probing Superexcited State Dynamics

Universal Probe of Reactive Scattering

XUV Probes of Surface Dynamics

(M. G. White)


Longer-term

Time-resolved photoelectron imaging


Nonlinear processes in the XUV

(L. DiMauro)

Installation of the first experiment - Dec, 2002




Monochromator





Ion Pair Imaging station

Ion Pair Imaging Spectroscopy



IPIS Technique: Excitation of a molecule in the VUV/XUV accesses ion-pair states that dissociate. If one of the products is structureless (such as F^-), then the kinetic energy release directly reflects the internal energy in the other product and the dynamics of the process.

Arthur Suits, BNL Chemistry and SBU



Ion Pair Dissociation of CH₃F


F⁻ Images

Photon
Energy:

13.52 eV

13.68 eV




13.95 eV



XUV
Polarization

DUV-FEL Upgrade - 100 uJ at 100 nm

FEL Scaling With NISUS & VISA



Summary

Tremendous progress was made in DUV-FEL in the last couple years.

We have developed a plan to upgrade DUV-FEL to 300 MeV and reach 100 nm for DUV-FEL.