

Hard X-ray Micro-imaging Beamlines at SPring-8

-Present status and some applications-

Yoshio Suzuki,
Microimaging/Analysis Group, SPring-8

There are four public beamlines used for micro-imaging in SPring-8, BL20XU, BL20B2, BL37XU, and BL47XU. Although, these beamlines are not exclusively used for micro-imaging, most of micro-imaging activities can be covered by combining these four beamlines. Details of these beamlines and some experimental results will be presented.

BL20XU is a medium-length (248 m) beamline with undulator light source [1,2]. Although, at the beginning, BL20XU was proposed as a medical imaging beamline, the present major fields in 20XU are R&D of optical systems and optical elements in hard X-ray region (8-110 keV), and coherent optics (holography, speckle, etc.), because its unique property is high-brilliance and high spatial-coherence by utilizing the undulator X-ray source and long beam-transport channel. An example of experimental results is shown in Fig. 1. Hard X-ray microbeam with 100 nm-resolution is attained with a FZP optics [2,3]. The optical devices and optics developed at BL20XU are transferred to other beamlines, and supplied to user experiments. Adding to them, a phase-contrast CT with Bonse-Hart type interferometer is now opened to public use at BL20XU, and 3-dimensional phase-contrast images can be taken with a spatial resolution of 10 μm [4].




Fig. 1. X-ray microbeam and scanning microscopy experiment at BL20XU. Outermost zone width of FZP is 100 nm, and focal length is 100 mm at 8 keV. (a) focused beam profile of 1st-order diffraction, (b) scanning microscopy image of test patterns, (c) focused beam profile of 3rd-order diffraction. Spatial resolution near to the diffraction-limit is achieved.

BL20B2 is a medium length (215 m) beamline with bending magnet light source[5]. This beamline is used in many fields of imaging techniques; micro-CT with spatial resolution of 10-100 μm , refraction-enhanced imaging, diffraction topography, micro-angiography, etc. The experimental setup for these application is routinely used by many user groups. Majority of the users in BL20B2 is micro-CT for material sciences, planetary sciences, and biomedical applications. The unique feature of BL20B2 is wide beam-cross-section. The beam width at the end station is 300 mm. A single-shot X-ray topography of 12-inches diameter Si wafer can be taken, as shown in Fig. 3. The high spatial-coherence and wide beam size is also appropriate for characterization of X-ray optical element with a large aperture, i.e. X-ray telescope. The 15 m-long experimental hutch at the end station is also significant for this purpose. Characterization (imaging properties and efficiency, etc.) of a hard X-ray telescope with 400 mm-aperture and 8 m-focal length is now being done at BL20B2.

BL37XU is a new public beamline constructed for X-ray fluorescence (XRF) analysis [6]. The BL37XU is designed to be a user-friendly beamline for spectro-microscopy, and high energy XRF experiment up to 75 keV can be done. Microbeam and scanning microscopy with spatial resolution of about 1 μm is routinely performed using the FZP objectives and/or total-reflection mirror optics in the energy region around 15 keV. The microbeam is now used by most of user groups at BL37XU.

BL47XU is named "R&D" beamline, and was designed as a pilot beamline in SPring-8. Now, about 40 % of total machine time of BL47XU is opened for general users, as a multipurpose beamline.

Major activity of public users is micro-CT and microbeam applications. Projection-type micro-CT with $1\ \mu\text{m}$ -resolution is now routinely used for material science and geophysics [7], and imaging-type CT with sub- μm resolution is also achieved using a FZP objective, as shown in Fig. 3 [8]. Micro-beam diffraction is also performed with $0.3\ \mu\text{m}$ -probe size, and now applied to material science [9].

We consider that nondestructive measurement by the combination of micro-CT and microbeam analysis (scanning microscopy, fluorescent X-ray spectroscopy, and micro-diffraction) will be an important tool for material sciences, planetary sciences, and biology.




Fig. 2. X-ray Topography of 300 mm-diameter Si wafer taken at BL20B2
X-ray energy: 21 keV




Fig. 3. X-ray microtomography using Fresnel zone plate objective at BL47XU
Sample: stony meteorite Allende (top), and diatom (bottom). X-ray energy: 8 keV,

1. Y. Suzuki, et al., Journal de Physique IV France 104 (2003) 35.
2. Y. Suzuki, et al., AIP Conference Proceedings 705 (2004) 344.
3. H. Takano, et al., Jpn. J. Appl. Phys. 43 (2003) L132.
4. A. Momose, et al., Journal de Physique IV France 104 (2003) 599.
5. S. Goto, et al., Nucl. Instrum. and Meth. A467-468 (2001) 682.
6. Y. Terada, et al., AIP Conference Proceedings 705 (2004) 376.
7. K. Uesugi, et al., Nucl. Instrum. and Meth. A467-468 (2001) 853.
8. A. Takeuchi, et al., Rev. Sci. Instrum. 73 (2002) 4246.
9. T. Iwata, et al., Macromol. Rapid. Commun. 24 (2004) 1100.