

Sub-psec Magnetic Domain Imaging in Nanoscale Structures Using Photoemission Electron Microscopy

David Keavney, Advanced Photon Source, Argonne National Laboratory

Advanced Photon Source

J.W. Freeland

R.A. Rosenberg

P. Ryan

R. Winarski

University of Utah

W. Casey Uhlig


Di Wu

Jing Shi




Use of the Advanced Photon Source is supported by the U.S. Department of Energy, Office of Science, under Contract No. W-31-109-ENG-38. Work at the University of Utah is supported by the National Science Foundation ECS under grant No. 0115164.

Time scales in magnetism




Applications of nanoscale magnets


Media




MRAM




Spin dynamics and time scales




Needed time resolutions:
< 200 ps




< 1 ps






< 1 ps



X-ray Photoemission Electron Microscopy at the APS




Chemical and Magnetic Imaging




- Simultaneous chemical + magnetic contrast
- High magnetic contrast
- Full field
- No interaction with sample
- 100-120 nm resolution
 - 20 nm target
- Domain imaging
- Buried layers (~5 nm)
- Finite size effects
- Self-assembled systems
- Ground states in nanoscale systems
- Interactions in particle arrays

Co nanodot arrays

e-beam lithographic lift-off process:
University of Utah








500x1000nm, 200x600nm








Magnetostatic interactions between Co dots

500 x 1000 nm dots






Magnetostatic interactions between Co dots

200 x 600 nm dots




Time resolved imaging


Pump-probe experiments using the APS ring timing structure



Magnetic linear dichroism in threshold photoemission





100 W Hg lamp
(4.9 eV)



0.37% asymmetry


Spatial Resolution at threshold

Chromatic aberrations limit the spatial resolution




Pd on Si

Spherical aberration limit ~ 20 nm




Magnetic imaging at 300 fs, 20-30 nm!

Beyond stroboscopic imaging



30s exp ($\sim 3 \times 10^{14}$ photons)



2-3% noise


4.9 eV @ 30 μJ = 3.8×10^{13} photons

200 μJ = 2.5×10^{14} photons

Enhancements:

Asymmetry (energy/polarization dependence?)

Flux (2mJ/pulse?)



Single-shot mode?

Summary

- Significant improvement in state of the art magnetic imaging possible with existing parameters
 - Magnetic nanostructures
 - Spin populations in semiconductors
- Possible single-shot imaging mode
 - Spin fluctuations