

# ***Using IDL and Python with EPICS***

*Mark Rivers, University of Chicago*

## ***Argonne National Laboratory***



*A U.S. Department of Energy  
Office of Science Laboratory  
Operated by The University of Chicago*



# Outline

---

- **Quick Overview of IDL**
- **ezca library**
- **Calling ezca from IDL**
- **IDL CA API**
- **IDL EPICS class libraries**
- **IDL applications**
- **Overview of Python**
- **Python class libraries**
- **Python applications**
- **Using EPICS from Visual Basic**

# Overview of IDL

---

- A high-level interpreted *programming language* with vector and array primitives - sort of a cross between BASIC and APL
- Modern programming language
  - Flow control
  - Data structures
  - Objects
- All operators and most functions work on scalar, vector or array data of any data type.
- Data *visualization tool*, advanced built-in graphics
  - 2-D plots
  - Contour plots
  - Surface plots
  - Shaded surfaces
  - Gray scale/color images
  - Isosurfaces
  - Volume rendering
- Multi-platform support
  - Unix: Sun, Hewlett Packard, Silicon Graphics, IBM
  - Linux
  - Microsoft Windows
  - Mac Darwin
- List price: ~\$3,000 on workstations, ~\$1,500 on PC, Mac

# Overview of IDL

---

- **Can call external C or other code**
- **Very fast for array operations, as fast as compiled languages**
- **GUI builder**
- **Multi-threaded**
- **Good vendor support**
- **IDL GUI applications can be run at no cost (IDL 6.0 and above)**
  - Must have license to use IDL command line

# Overview of IDL

---

## Data Structures

- **A variable in IDL has both a structure and a data type associated with it. Both of these are dynamic, i.e. they can be changed via an assignment statement at any time.**
- **Data types**
  - Byte (b bit, unsigned)
  - Integer (16 bit, signed)
  - Long (32 bit, signed)
  - Float (32 bit floating point)
  - Double (64 bit floating point)
  - Complex (pair of 32 bit floats)
  - Double complex (pair of 64 bit floats)
  - String (0 to 64k characters)
- **Data Structures**
  - Scalar
  - Vector
  - Array - up to 7 dimensions
  - Structure - composed of other elements, like C
- **The sizes of arrays are limited only by the available virtual memory.**

# Overview of IDL

---

## Assignment Statements

`A = B + 1`

- **A has the same structure as B, with a data type equal to that of the most precise operand in the expression on the right hand side. In this case it could be any type except string.**
- **If B is a vector or array then 1 is added to each element.**

`A = 0` ; A is a 16 bit integer

`A = A * 0.5` ; A is now a 32 bit float

`B = A(*,3)` ; B is equal to the 4th row of A

`A(*,3) = 0` ; Set all elements in 4th row of A equal to 0

## Syntax

- **Examples:**

`image = fltarr(512, 512)` ; zero filled array

`b = image(0:127, 0:127)` ; b is 128x128 array

`image(*,100) = findgen(512)` ; replace row 100

`plot, image(*,120)` ; plot row 121

; Display the power spectrum as an image

`tvscrl, alog(abs(fft(image, 1)))`


# IDL Examples

```
IDL> a = sin(findgen(100)/99. * 2 * !pi)
```

```
IDL> help, a
```


```
A                FLOAT      = Array[100]
```

```
IDL> plot, a
```



# IDL Examples

```
IDL> a = shift(alog(abs(fft(dist(256),1))),128,128)  
IDL> isurface, a
```





# ezca and EzcaScan

- **EPICS extensions for “Easy Channel Access”**

- Don't need to handle chids, just use PV name strings – hash table
- Synchronous APIs - applications don't have to handle callbacks

- **Ezca (partial list)**

```
- epicsShareFunc int epicsShareAPI ezcaGet(char *pvname, char ezcatype,  
-     int nelem, void *data_buff);  
- epicsShareFunc int epicsShareAPI ezcaPut(char *pvname, char ezcatype,  
-     int nelem, void *data_buff);  
- epicsShareFunc int epicsShareAPI ezcaPutOldCa(char *pvname, char ezcatype,  
-     int nelem, void *data_buff);  
- epicsShareFunc int epicsShareAPI ezcaNewMonitorValue(char *pvname,  
-     char ezcatype); /* returns TRUE/FALSE */  
- epicsShareFunc int epicsShareAPI ezcaSetTimeout(float sec);  
- epicsShareFunc float epicsShareAPI ezcaGetTimeout(void);  
- epicsShareFunc int epicsShareAPI ezcaSetRetryCount(int retry);  
- epicsShareFunc int epicsShareAPI ezcaGetRetryCount(void);  
- epicsShareFunc int epicsShareAPI ezcaPvToChid(char *pvname, chid **cid);  
- epicsShareFunc int epicsShareAPI ezcaSetMonitor(char *pvname, char ezcatype);  
- epicsShareFunc int epicsShareAPI ezcaClearMonitor(char *pvname, char ezcatype);  
- epicsShareFunc int epicsShareAPI ezcaStartGroup(void);  
- epicsShareFunc int epicsShareAPI ezcaEndGroup(void);  
- epicsShareFunc int epicsShareAPI ezcaGetControlLimits(char *pvname,  
-     double *low, double *high);  
- epicsShareFunc int epicsShareAPI ezcaGetGraphicLimits(char *pvname,  
-     double *low, double *high);  
- epicsShareFunc int epicsShareAPI ezcaGetNelem(char *pvname, int *nelem);  
- epicsShareFunc int epicsShareAPI ezcaGetPrecision(char *pvname,  
-     short *precision);  
- epicsShareFunc int epicsShareAPI ezcaGetStatus(char *pvname,  
-     TS_STAMP *timestamp, short *status, short *severity);  
- epicsShareFunc int epicsShareAPI ezcaGetUnits(char *pvname,  
-     char *units); /* units must be at least EZCA_UNITS_SIZE large */
```

# *ezca and EzcaScan*

---

- **EzcaScan (partial list)**

- `epicsShareFunc int epicsShareAPI Ezca_getArray(int noName, char **pvName, int type, int nodata, void *value);`
- `epicsShareFunc int epicsShareAPI Ezca_getArrayEvent(int noName, char **pvName, int type, int nodata, void *value);`
- `epicsShareFunc int epicsShareAPI Ezca_putArray(int noName, char **pvName, int type, int nodata, void *value);`
- `epicsShareFunc int epicsShareAPI Ezca_putArrayEvent(int noName, char **pvName, int type, int nodata, void *value);`

# ***ezca and IDL***

---

- **IDL can call “shareable libraries”, e.g. .so files on Unix, .dll files on Windows**
- **The argument passing convention is fixed, it is not compatible with ezca.dll directly**
- **Need a thin glue layer between IDL and ezca/EzcaScan**
- **ezcaIDL is the glue layer. Mostly just changes calling conventions. Provides a few functions that ezca and EzcaScan do not. Use the ezcaPvToChid() function.**
  - **ezcaIDLGetEnumStrings()**
  - **ezcaIDLGetCountAndType()**

# ***IDL Channel Access API***

---

## **Routines which return information about process variables**

```
Status = caGet(pvname, value, /string, maximum=max)
Status = caGetControlLimits(pvname, low, high)
Status = caGetGraphicLimits(pvname, low, high)
Status = caGetPrecision(pvname, precision)
Status = caGetStatus(pvname, timestamp, status, severity)
Status = caGetUnits(pvname, units)
Status = caGetEnumStrings(pvname, strings)
Status = caGetCountAndType(pvname, count, type)
```

## **Routines which write new values to process variables**

```
Status = caPut(pvname, value, wait=wait)
```

## **Routines which control channel access timeouts**

```
Timeout = caGetTimeout()
caSetTimeout, timeout
RetryCount = caGetRetryCount()
caSetRetryCount, retrycount
```

# ***IDL Channel Access API***

---

## **Routines which control synchronous groups**

`caStartGroup`

`stat = caEndGroup(status)`

## **Routines which control channel access monitors**

`Status = caSetMonitor(pvname)`

`Status = caClearMonitor(pvname)`

`State = caCheckMonitor(pvname)`

## **Routines which control debugging and error messages**

`caDebug, state`

`caTrace, state`

`caError, err_string, /ON, /OFF, /PRINT, prefix=prefix`



## **Documentation at**

<http://cars.uchicago.edu/software/idl/ezcaIDLGuide.html>

<http://cars.uchicago.edu/software/idl/ezcaIDLRef.html>

# IDL EPICS Examples

```
IDL> status = caget('13LAB:m1.VAL', position)
IDL> help, status, position
STATUS          LONG          =          0
POSITION        DOUBLE        =      517.19305
IDL> status = caget('13LAB:quadEM:mcal', spectrum)
IDL> plot, spectrum
IDL> help, status, spectrum
STATUS          LONG          =          0
SPECTRUM        LONG          = Array[2048]
IDL> plot, spectrum[0:500]
IDL> fft_data = alog(abs(fft(spectrum,1)))
IDL> plot, fft_data, xrange=[0,1023]
```



# ***IDL EPICS Examples***

---

## **Move a motor**

```
DL> status = caput('13LAB:m8.VAL', 10000)
IDL> status = caget('13LAB:m8.RBV', pos)
IDL> print, pos
      215.52734
IDL> status = caget('13LAB:m8.RBV', pos)
IDL> print, pos
      835.64453
IDL> status = caget('13LAB:m8.RBV', pos)
IDL> print, pos
     1795.6055
```

## **Ezca timeout values are important!**

```
IDL> print, cagettimeout()
      0.0200000
IDL> t0=systime(1)&for i=1,100 do t=caget('13LAB:m1', v)&print, systime(1)-t0
      2.9898720
IDL> casetimeout, .001
IDL> t0=systime(1)&for i=1,100 do t=caget('13LAB:m1', v)&print, systime(1)-t0
      0.21649790
```

# IDL EPICS Examples

---

## Using monitors

Monitored channels read the cached values on `caget()`

Can check whether a monitor has happened (a Channel Access value callback)

```
IDL> status = caSetMonitor('13LAB:m8.DMOV')
IDL> state = caCheckMonitor('13LAB:m8.DMOV')
IDL> help, state
STATE          LONG          =          1
IDL> status = caget('13LAB:m8.DMOV', done)
IDL> help, done
DONE           INT           =          1
IDL> state = caCheckMonitor('13LAB:m8.DMOV')
IDL> help, state
STATE          LONG          =          0
IDL> status = caput('13LAB:m8.VAL', 0)
IDL> state = caCheckMonitor('13LAB:m8.DMOV')
IDL> help, state
STATE          LONG          =          1
IDL> status = caget('13LAB:m8.DMOV', done)
IDL> help, state
STATE          LONG          =          1
IDL> help, done
DONE           INT           =          0
```

**Monitors are useful for seeing that a PV changed state, even if its value is the same because one “missed” the transition. For example, PV goes 0->1->0. IDL polling might miss the one state, but checking a monitor would let one know that it happened.**



# ***IDL EPICS Class Libraries***

---

- **IDL object classes that hide the underlying EPICS process variables**
- **IDL objects treat all data as private, only accessible through methods.**
- **Provide an object-oriented interface to common beamline objects (motors, scalers, mcas, scans)**
  - `epics_motor`
  - `epics_scaler`
  - `epics_mca` (inherits device-independent mca class)
  - `epics_med` (multi-element detector)
  - `epics_sscan`
- **Example of `epics_motor`**

```
IDL> motor = obj_new('EPICS_MOTOR', '13LAB:m8')
```


```
IDL> motor->move, 1000.      ; Move to absolute position 10.
```

```
IDL> motor->wait              ; Wait for it to get there
```

# ***IDL EPICS Class Libraries***

- **Example of epics\_mca**

```
IDL> mca = obj_new('epics_mca', '13LAB:aim_adc1')
IDL> mca->erase
IDL> mca->acquire_on
IDL> data = mca->get_data()
IDL> plot, data
```



- **Example of epics\_scaler**

```
IDL> scaler = obj_new('epics_scaler', '13LAB:scaler1')
IDL> scaler->start, 10.      ; Count for 10 seconds
IDL> scaler->wait            ; Wait for it to get done
IDL> counts = scaler->read(); Read the counts on each channel
IDL> print, counts
100000000      0      0      0      0      0      0      0
```

# ***IDL EPICS Class Libraries***

---


- **epics\_sscan**
- **Designed to do the following:**
  - Provide an object-oriented interface to standard EPICS scans, enabling user written software to easily access scan header information and data.
  - Provide an easy way to read MDA files written by the saveData function in synApps.
  - Provide an easy way to get scan data into the IDL iTools system. iTools provide powerful interfaces for visualizing data, zooming in, adding annotation, and producing publication quality plots.
  - Provide a way to convert binary scan files (e.g. MDA) into ASCII
  - Does not currently communicate with the IOC for real-time data, but this is planned for the future

# IDL EPICS Class Libraries

## Example: Simple 1D epics\_sscan

```
IDL> s = read_mda('13IDC_0027.mda') ; Read the data
```

```
IDL> s->display ; Display the first detector
```




# IDL EPICS Class Libraries

## Example: 2-D epics\_sscan

IDL> s=read\_mda('2idd\_0087.mda') ; Read the 2-D dataset

IDL> s->display, /all, /grid ; Display all of the images in a grid




# IDL EPICS Class Libraries

## Example: 2-D epics\_sscan

; Plot a profile of column 20 (X=20) in detector 15.

```
IDL> s->display, detector=15, xrange=20
```



# IDL EPICS Class Libraries

## Documentation: Reference manual for each class library

### EPICS\_MCA Class

This page was created by the IDL library routine `mk_html_help`. For more information on this routine, refer to the IDL Online Help Navigator or type:

?mk\_html\_help at the IDL command line prompt.

Last modified: Sat Jul 14 10:16:05 2001.

### List of Routines

- [EPICS\\_MCA::ACQUIRE\\_OFF](#)
- [EPICS\\_MCA::ACQUIRE\\_ON](#)
- [EPICS\\_MCA::ACQUIRE\\_WAIT](#)
- [EPICS\\_MCA::ADD\\_ROI](#)
- [EPICS\\_MCA::DEL\\_ROI](#)
- [EPICS\\_MCA::ERASE](#)
- [EPICS\\_MCA::GET\\_ACQUIRE\\_STATUS](#)
- [EPICS\\_MCA::GET\\_CALIBRATION](#)
- [EPICS\\_MCA::GET\\_DATA](#)
- [EPICS\\_MCA::GET\\_ELAPSED](#)
- [EPICS\\_MCA::GET\\_PRESETS](#)
- [EPICS\\_MCA::GET\\_ROIS](#)
- [EPICS\\_MCA::GET\\_ROI\\_COUNTS](#)
- [EPICS\\_MCA::GET\\_SEQUENCE](#)
- [EPICS\\_MCA::INIT](#)
- [EPICS\\_MCA::SET\\_CALIBRATION](#)
- [EPICS\\_MCA::SET\\_DATA](#)
- [EPICS\\_MCA::SET\\_PRESETS](#)
- [EPICS\\_MCA::SET\\_ROIS](#)
- [EPICS\\_MCA::SET\\_SEQUENCE](#)
- [EPICS\\_MCA::SPECTRA\\_SCAN](#)
- [EPICS\\_MCA::WRITE\\_FILE](#)
- [EPICS\\_MCA\\_DEFINE](#)
- [RELEASE NOTES](#)

Documentation at:

<http://cars9.uchicago.edu/software/idl/>

#### EPICS\_MCA::GET\_DATA

[\[Previous Routine\]](#) [\[Next Routine\]](#) [\[List of Routines\]](#)

##### NAME:

EPICS\_MCA::GET\_DATA

##### PURPOSE:

This function returns the data from the MCA.

##### CATEGORY:

EPICS device class library.

##### CALLING SEQUENCE:

Result = epics\_mca->GET\_DATA()

##### KEYWORD PARAMETERS:

###### CHECK\_NEW:

A flag which indicates that this routine should only return the data if it has changed.

##### OPTIONAL OUTPUTS:

###### NEW\_FLAG:

If CHECK\_FLAG is set, then NEW\_FLAG will be 1 if the function is returning new data, 0 if the function is not returning new data. If CHECK\_FLAG is set and NEW\_FLAG is 0 then the function returns -1.

##### PROCEDURE:

This function reads the data from the hardware using the EPICS MCA record, and then invokes MCA::GET\_DATA

##### ADDITIONAL INFORMATION:

See [MCA::GET\\_DATA\(\)](#).



##### MODIFICATION HISTORY:

Written by: Mark Rivers, October 1, 1997  
Nov. 14, 1997 Mark Rivers. Changed routine to eliminate setting rec.READ back to 0, since record support does this automatically and it was causing record to process again.  
19-Sep-1998 MLR Added /WAIT to caput, since default is not to wait for callback now.  
17-Mar-1999 MLR Removed /WAIT from caput, to be compatible with version 4.3 and later of the MCA record, which does not fire forward links until acquisition is complete.  
28-Mar-1999 MLR Changed routine so it no longer pokes READ field. This assumes that someone else (typically a database) is periodically poking the READ field. The object initialization code now sets a monitor on the VAL field. Added New\_flag output and CHECK\_NEW keyword.

(See `epics_mca_define.pro`)

# IDL EPICS Applications




- **mcaDisplay**
  - Full-featured program for displaying, controlling EPICS multi-channel analysers, including peak fitting
  - Uses epics\_mca class library, and exports mca\_display class, so it can be controlled by other IDL applications





# IDL EPICS Applications

## Data catcher and data viewer (Ben-Chin Cha)



# Using EPICS from Visual Basic

---

- **ezca.dll** can be called directly from Visual Basic on Windows
- **ezca.bas** provides the interface

```
Public Const ezcaByte As Byte = 0
Public Const ezcaString As Byte = 1
Public Const ezcaShort As Byte = 2
Public Const ezcaLong As Byte = 3
Public Const ezcaFloat As Byte = 4
Public Const ezcaDouble As Byte = 5
```

```
Public Declare Function ezcaGet Lib "ezca.dll" _
    (ByVal pvname As String, _
    ByVal ezcatype As Byte, _
    ByVal nelem As Long, _
    ByRef data As Any) As Long
```

```
Public Declare Function ezcaPut Lib "ezca.dll" Alias "ezcaPutOldCa" _
    (ByVal pvname As String, _
    ByVal ezcatype As Byte, _
    ByVal nelem As Long, _
    ByRef data As Any) As Long
```

```
Public Declare Function ezcaPutString Lib "ezca.dll" Alias "ezcaPutOldCa" _
    (ByVal pvname As String, _
    ByVal ezcatype As Byte, _
    ByVal nelem As Long, _
    ByVal data As Any) As Long
```

```
Public Declare Function ezcaPutCallback Lib "ezca.dll" Alias "ezcaPut" _
    (ByVal pvname As String, _
    ByVal ezcatype As Byte, _
    ByVal nelem As Long, _
    ByRef data As Any) As Long
```

# Using EPICS from Visual Basic

- Example: tomography data collection. VB used because it can easily control Roper's WinView program for the CCD detector

