




# Agenda

- Discuss group objectives  
Create a forum/discussion group to share information on
  - Hardware
    - Beam line characteristics
    - CCD cameras/Scintillator screens
    - Sample mounting
    - Stages
  - Software
    - Experiment control software
    - Reconstruction code
    - Data: format, archival, distribution ...
  - Problems
  - User sharing
- Define and document available tomography resources at the APS
  - Define: Presentations
  - Document: Connect web pages ?
- Future meeting schedules
- Presentation: 2BM setup


# 2BM-A










# Tomography System




# Tomography Setup



- Scintillator

- CdWO<sub>4</sub> 10 x 10 x 0.5 mm polished
- CdWO<sub>4</sub> 10 x 10 x 0.3 mm cleaved
- CdWO<sub>4</sub> 25 x 25 x 0.5 mm
- YAG:Ce film 1  $\mu\text{m}$  on undoped YAG 10x10 mm x 100  $\mu\text{m}$  substrate
- YAG:Ce film 5  $\mu\text{m}$  on undoped YAG 10x10 mm x 300  $\mu\text{m}$  substrate





# Tomography Experiment Control Software Functions





Automatically collect and reconstruct a tomography data set. This includes:

- *Setup (beam line specific parameters, sample, GU proposal, user info)*
- *Data Acquisition (control of CCD camera, sample stages, robot, beamline shutter)*
- *Preprocessing*
- *Reconstruction (cluster configuration)*
- *Postprocessing*

Generate an electronic log for the experiment



# Setup



Manufacturer	Model	Comment
Zeiss	1.25x	5.08 um/pixel
Zeiss	4x	1.66 um/pixel
Zeiss	5x	1.33 um/pixel
Zeiss	10x	0.66 um/pixel
Zeiss	20x	0.33 um/pixel
Zeiss	Tube Lens	Standard 1x AXIOPLAN mic...





  

Add			
Order	Manufacturer	Model	Comment
1	Zeiss	Tube Lens	Standard 1x AXIOPLAN...
2	Zeiss	5x	1.33 um/pixel


  

Remove			
OK	Cancel	OK	Cancel



# Setup



# Data Acquisition



# Reconstruction



# Generate an electronic log for the experiment


The screenshot shows a software interface for managing an HDF5 file named 'test.hdf'. The left pane, titled 'Groups', displays the hierarchical structure of the file. The root group 'test.hdf' contains 'facility' and 'user' groups. Below 'test.hdf' is the 'experiment' group, which is expanded to show its subgroups: 'setup', 'storage\_rir', 'beamline', 'sample', 'motors', 'scintillator', 'objectives', and 'detector'. The 'detector' group is highlighted with a yellow selection bar. The right pane, titled 'Fields', lists the metadata fields and their values. The table has columns for File Index, Data Type, Rank, Size, and Value.

File Index	Data Type	Rank	Size	Value
prefix	char	1	4	2bmb
device	char	1	5	csccd
roi	char	1	2	01
manufacturer	char	1	16	Roper Scientific
model	char	1	7	Unknown
bit_depth	float32	1	1	16.0
exposure_time	float32	1	1	1.01
bin_x	int32	1	1	1
bin_y	int32	1	1	1
size_x	int32	1	1	1024
size_y	int32	1	1	1024
orientation	float32	1	1	90.0
magnification	float32	1	1	1.0
comment	char	1	23	This is a CCD ...

## C:\TomoMPI

File View Analyse Tools Help Special

Cross Correlation Shift (X, Y, Xshift)	0	0	0
Start	0	Weight X Shift	0
End	0	Weight Tilt	0
		Manual X Shift	0
		Manual Tilt Angle	0

 Rotate by 90 deg

**Experiment Name** S2\_Load0\_Y1.hdf  
**File Size** 524288  
**Loaded Bytes** 262144  
**Picture Size** ( 512 512 )  
**Cursor Position** ( 507 309 )  
**Pixel Value** 0.00194415438454598  
**Mean, StdDev, Var** ( 0.2999650.3606970.2051714 )

Data File 1  
Normalized Data File 3  
Sinogram File 1  
Reconstruction 1

Showing Normalized Data D:\Load0\_Y1\normalized\norm\_S2\_Load0\_Y1\_00003.hdf at 0 deg

### ExperimentOpenFile

Experiment Selection  
File Name D:\test.hdf

**Proposal Number:** GU-323  BIN

User  
Name Stock, Stuart  
Affiliation Northwestern Univ. Medical School  
Address 303 E. Chicago Ave., Chicago, IL  
Phone (708) 756 8324  
E-mail s-stock@northwestern.edu

Sample  
Name test2048  
Description test data set  
Confidential NO  
Comment Somthing wicked this way comes...  
Data Collected On Sun Nov 24 14:15:31 CST 2002

Facility  
Name APS Sector 2 Beamline BM Station B

Setup

Beamline

Beam Mode	Mono
Filter A	None
Filter B	None
A Slits X/Y (Center/Width)	0 6 0 1.6
Mirror Coating	Cr
Monocromator	Kohzu DCM 15
B Slits X/Y (Center/Width)	0 1.2 32 1





CCD Camera





Name	2bmb csccd 01 Roper	CoolSnap
Bit Depth	16	
Exposure Time (s)	6	
Bin (X, Y)	1 1	
Size (X, Y)	1024 1024	
Objective 1	Zeiss	Tube Lens
Objective 2	Zeiss	5x
Scintillator	YAG	YAG
Magnification	1	
Comment	Somthing wicked this way comes...	


Acquisition Parameters




Angle (Start, End, interval)	0 180 0.25
White/Dark Field Frequency	20
Total Files	797

Experiment File Created On Thu Jan 16 15:38:23 CST 2003









# What we are developing:

- **Phase-Contrast Tomography Techniques**
  - Reconstruction algorithms
  - Local Tomography Problem
- **Sample Changer integration in the beam line software**

# What we don't have:

- **Custom 3D rendering tools**
- **Post reconstruction data analysis**
- ...