

# ALS/APS Collaboration on 480 x 480 Fast CCD Detector

## The APS story


J. Weizeorick  
Beamline Technical Support Group  
Advanced Photon Source, X-ray Science Division  
Argonne National Laboratory

**2009 Advanced Light Source Users' Meeting**  
October 15-17, 2009

# Outline

- Beamline Technical Support Group at APS  
(Slides from P. Fernandez)
- APS 480 x 480 Fast CCD Story
  - Milestones
  - Scientific Data (Slides from Dr. Larry Lurio from NIU)
  - Improvements
  - Detector Specifications
  - Future Tasks and Goals
- Future of Fast CCD Detectors at ANL
  - 1Kx1K Frame Transfer CCD

# Advanced Photon Source Organization Chart



# The XSD Beamline Technical Support Group

- Two main responsibilities
  - The APS Detector Pool
  - X-ray detector development





# The APS Detector Pool

The **goals** of the Detector Pool are:

- To provide a range of detectors suitable for the research interests of the APS community.
- To provide expertise and guidance on x-ray detector capabilities and use.
- To identify and assess upcoming commercial technologies.
- Operate visible light and x-ray detector test stands, available to APS scientists.



APS Detector Pool – Contact: A. Miceli




# XSD-BTS X-ray Detector R&D Activities

## *ALS/APS Collaboration on 480 x 480 fast CCD detectors*

*P. Denes, H. Padmore, ALS; J. Weizeorick, APS/XSD-BTS*

- Produce two Fast CCD x-ray detectors using 480x480 pixels CCD chip, and six fCRIC readout chips
- At the APS, strong interest from the X-ray Photon Correlation Spectroscopy (XPCS) community
- BTS Personnel that worked on project
  - Jonathan Baldwin
  - Kevin Beyer
  - Lisa Gades
  - John Lee
  - Troy Lutes
  - Tim Madden
  - Antonino Miceli
  - Steve Ross
  - Chris Piatek




# APS 480 x 480 Fast CCD Story - 2009 Milestones






# APS 480 x 480 Fast CCD Story - 2009 Milestones





- January 11 – Peter and Dionisio arrived to install subassembly
- January 12 - First light with old substrate
- January 13 - First light with new substrate
- January 14 - First x-ray with x-ray tube
- January to July
  - Make Improvements
  - Characterize detector
- July 15, 2009 1<sup>st</sup> beam at 8id
- Nov 11, 2009 2<sup>nd</sup> beam at 8id






# APS 480 x 480 Fast CCD Story - Improvements


- January to July 2009 - Improvements and characterizing the detector
  - Improvements
    - Improved offset oscillation seen in a series of images
    - Fixed left to right intensity variation in 2x mode
  - Characterize Detector



# APS 480 x 480 Fast CCD Story - Science Data




ALS Fe fluorescence with 30 V Bias  
in 8x mode




APS Fe fluorescence with 30 V Bias  
in 8x mode

# APS 480 x 480 Fast CCD Story - Science Data



Histogram of ALS fluorescence data looking at bottom ½ of CCD





Histogram of APS fluorescence data looking at bottom ½ of CCD

S4 histogram algorithm - Sums the 4 most energetic pixels in a 5x5 box around a seed pixel. Requires that 70% of the energy in the box be in the 4 most energetic pixels.

# APS 480 x 480 Fast CCD Story - Science Data

The lens contains protein mixtures, that vary in concentration and composition with position.






# APS 480 x 480 Fast CCD Story - Science Data

## Scientific Motivation for Eye-Lens Protein Studies

In a healthy eye-lens the proteins are in a mixed state which yields the optimized index of refraction profile. Changes in the phase of the eye-lens proteins can lead to various forms of eye disease.

Cold cataract is due to reversible liquid-liquid phase separation in young, mammalian eye lenses.




A liquid-glass transition in eye-lens protein mixtures may be related to the stiffening of the eye lens seen in presbyopia, a loss of ability to focus the eye on nearby objects

# APS 480 x 480 Fast CCD Story - Science Data


X-ray Photon Correlation Spectroscopy (XPCS) can be used to measure both the structure and dynamics in mixtures of eye-lens proteins.

In the XPCS technique the sample is illuminated with a coherent x-ray beam which leads to the creation of a speckle pattern superimposed on the regular x-ray scattering pattern.

Dynamics within the sample cause the speckle pattern to fluctuate in time. Time autocorrelation of this fluctuating speckle pattern can be analyzed to yield information about the dynamics within the sample, such as localized diffusion coefficients.




Example of a speckle pattern produced using a coherent x-ray beam. (In this case a static aerogel sample)



Time correlation function from a concentrated alpha crystallin suspension. For this sample the correlation function remained at 1, which means that the dynamics were too fast to be successfully resolved with the available camera.


# APS 480 x 480 Fast CCD Story - Science Data

## Small-Angle X-ray Scattering Results



# Detector Improvements - Intensity oscillation

Images taken at APS X-ray tube in building 401 – 20090115\Tiff\BeamFrof\_400ms\*.tif




- Intensity Oscillation over whole image
- Bright rectangle is beam framed by slits
- Darker middle caused by two layers of electrical tape used to block out light
- Glow areas caused by hot spots on CCD




# Detector Improvements - Intensity oscillation

- Z Projection through a single pixel over image set - 20090219-FCCD\Dark\_Osc\_5msec\_exp-30deg\*.tif files
- Pattern repeating every 8 images
- Pixels 48,255 or pixels in column 1 have the biggest oscillation around 50 ADU
- Clock Phase problem between 125MHz(8nsec) and 40MHz(25nsec)


For Pixel (45,255)




For Pixel (46,255)



For Pixel (47,255)



For Pixel (48,255)



# Detector Improvements - Intensity oscillation

## - Original Configuration

**FCCD Interface Module**

**Detector ON-Line**

**CCD ON**

**CCD OFF**

**Clock Status**

**ON**

**Trigger Exposure(s)**

Exposure Mode: Focus

Exposure Time (CCD integration time): 100 msec

Num of Exps: 1

Exposure Cycle Time: 500 msec

Delay between T -> E: 0 usec

Shutter Time: 1 msec

Delay between E -> S: 0 usec

**Serial Port Status**

Port: COM2 Open

Baud: 115200 Close

**State Parameters**

	S0	S1	S2	S3	S4	S5	S6	S7
25nsec_per_tick	50	06	06	03	05	03	03	03
ticks_per_state	08	10	01	10	08	10	10	10
passes_per_state	01	0c	01	0c	01	07	05	07
next_state	01	00	05	04	07	06	04	00
loops_per_state	00	00	00	00	f9	00	00	00
loop_state	02	02	00	00	03	00	00	00

S0 - Vertical Shift during idle  
S1 - Horizontal Shift during idle  
S2 - Exposure - after S0 or S1 finishes  
S3 - Horizontal Read, convert and data saved  
S4 - Vertical Shift during read  
S5 - Horizontal Read Clocks with convert, data not saved  
S6 - Horizontal Read Clocks with convert, and data saved  
S7 - Horizontal Read Clocks with convert and data saved

**Clock Waveforms**

Hclk1-idle: 1111111000111111

Hclk2-idle: 0000001110000000

Hclk3-idle: 1111100011111111

Sw-idle: 0000000000000100

Rg-idle: 1111100000111111

Vclk1-idleRd: 1111000111111111

Vclk2-idleRd: 0001110000000000

Vclk3-idleRd: 1100011111111111

Tg-idleRd: 1000001111111111

Hclk1-Rd: 1111111000111111

Hclk2-Rd: 0000001110000000

Hclk3-Rd: 1111100011111111

Sw-Rd: 0000000000000001

Rg-Rd: 1111100000111111

Conv-Rd: 0000000001000000

**Voltages**

	+Clk Voltage 0V->+10V	-Clk Voltage -10V->0V		
V(1->3)	6.00	-2.00	OTG (0->5V)	3.00
TgClk	6.00	-2.00	VDDRST (0->-15V)	-15.000
H(1->3)	8.00	-3.00	VDDOUT (0->-25V)	-25
SwClk	5.00	-5.00	NGD (-10V->+10V)	0
RgClk	0.00	-6.00	NCON (-10V->+10V)	0
			GAURD (-10V->+10V)	0

Command Out: 41

Hex Command: 41

Reset ICRCs

Reset Data Modules Clock Sync

DataModule1: f0 00

DataModule2: f0 10

Powerup Mask = 00 00 00 00

DM1 Readout Mode: Descramble Right to Left

DM2 Readout Mode: Descramble Right to Left

Select Command File: Send Commands from File

C:\Documents and Settings\weizeorick\My Documents\Visual Studio 2008

ReadClkPtns: HELP: PURGE TEXT

MESSAGE BOX: Exposure time sent out 100

# Detector Improvements - Intensity oscillation

- Fixed on 3-02-2009 changes S2-25nsec\_per\_tick from 6 to 7

**FCCD Interface Module**

**Detector ON-Line**

**Trigger Exposure(s)**

**CCD ON**

**CCD OFF**

**Clock Status**

**ON**

Exposure Mode: Single

Exposure Time (CCD integration time): 900 msec

Num of Exps: 1

Exposure Cycle Time: 1000 msec

Delay between T -> E: 0 usec

Shutter Time: 200 msec

Delay between E -> S: 0 usec

**Serial Port Status**

Port: COM2 Open

Baud: 115200 Close

**Clock Waveforms**

Hclk1-idle: 1111111000111111

Hclk2-idle: 0000001110000000

Hclk3-idle: 1111100011111111

Sw-idle: 0000000000000100

Rg-idle: 1111100000111111

Vclk1-idleRd: 1111000111111111

Vclk2-idleRd: 0001110000000000

Vclk3-idleRd: 1100011111111111

Tg-idleRd: 1000001111111111

Hclk1-Rd: 1111111000111111

Hclk2-Rd: 0000001110000000

Hclk3-Rd: 1111100011111111

Sw-Rd: 0000000000000001

Rg-Rd: 1111100000111111

Conv-Rd: 0000000001000000

**State Parameters**

	S0	S1	S2	S3	S4	S5	S6	S7
25nsec_per_tick	50	06	07	03	05	03	03	03
ticks_per_state	08	10	01	10	08	10	10	10
passes_per_state	01	0c	01	0c	01	07	05	07
next_state	01	00	05	04	07	06	04	00
loops_per_state	00	00	00	00	f9	00	00	00
loop_state	02	02	00	00	03	00	00	00

S0 - Vertical Shift during idle  
S1 - Horizontal Shift during idle  
S2 - Exposure - after S0 or S1 finishes  
S3 - Horizontal Read, convert and data saved  
S4 - Vertical Shift during read  
S5 - Horizontal Read Clocks with convert, data not saved  
S6 - Horizontal Read Clocks with convert, and data saved  
S7 - Horizontal Read Clocks with convert and data saved

**Voltages**

+Clk Voltage -Clk Voltage	OTG (0->5V)	VDDRST (0->-15V)	VDDOUT (0->-25V)	NGD (-10V->+10V)	NCON (-10V->+10V)	GAURD (-10V->+10V)
0V->+10V -10V->0V	BIAS (0->99V)					
V(1->3): 6.00 -2.00	80.00	3.00	-15.000	0	0	0
TgClk: 6.00 -2.00						
H(1->3): 8.00 -3.00						
SwClk: 5.00 -5.00						
RgClk: 0.00 -6.00						

Command Out: Hex Command 41

Reset ICRCs

Reset Data Modules Clock Sync

DataModule1 DataModule2

Update Masks: f0 00 f0 10

Powerup Mask = 00 00 00 00

DM1 Readout Mode: Descramble Right to Left

DM2 Readout Mode: Descramble Right to Left

Select Command File Send Commands from File

C:\Documents and Settings\weizeorick\My Documents\Visual Studio 2008


ReadClkPtns HELP PURGE TEXT

MESSAGE BOX: Exposure time sent out 900


# Detector Improvements - Intensity oscillation

- Z Projection through a single pixel over image set
  - 20090310-FCCD\DarkLightsOn-30C-10msec\DarkLightsOn\_10msecE\_1sE-30C\*.tif files
- Cannot see a pattern
- Column 1 went from around a 50ADU oscillation to around a 10 ADU oscillation


For Pixel (45,255)




For Pixel (46,255)



For Pixel (47,255)




For Pixel (48,255)






# Detector Improvements - Right to Left Intensity variation

- Over scan removed between every 10 columns




# Detector Improvements - Right to Left Intensity variation

- Photon Transfer Curve




Fast CCD Detector with integrating sphere

Movie of photon transfer curve Data Set



# Detector Improvements - Right to Left Intensity variation


- Photon Transfer Curve




From: Scientific Charge-Coupled Devices  
By James R. Janesick

# Detector Improvements - Right to Left Intensity variation

- Photon Transfer Curves Variation from Column 1 to Column 10
- Spent long time with analysis
- Eventually moved on to test out 8x mode
  - 8x mode didn't have the right to left intensity variation!



Column 1




Column 10




# Detector Improvements - Right to Left Intensity variation

- Made configuration of 2x Mode like 8x mode on 4-07-09
  - Changed “25nsec\_per\_tick” in S3, S5, and S6 from 03 to 04
    - S3 = horizontal read, convert and save data
    - S5 = horizontal read, convert and don't save data (flush fCRIC pipe)
    - S6 = horizontal read, convert and save data (flush finished)
  - S5 “passes\_per\_state” changed from 7 to 6
  - S6 “passes\_per\_state” changed from 5 to 6




# Detector Improvements - New 2x Configuration



Dark image with new  
2x configuration


# Detector Improvements - New 2x Configuration



Dark image with new  
2x configuration


Same image  
-Over Scan Removed  
-B&C adjusted

# Detector Improvements - New 2x Configuration



Dark image with new  
2x configuration


Same image  
-Over Scan Removed  
-B&C adjusted




Horizontal cut of image

# Detector Improvements - New 2x Configuration

- Photon Transfer Curves Same for Column 1 to Column 10



Column 1



Column 10



# APS 480 x 480 Fast CCD Story - Detector Specifications


- 480 x 480 Pixels
- 30 x 30  $\mu\text{m}$  pixels size
- 125 Frames Per Second
- 5.1 msec readout

14.4 mm x 14.4 mm CCD area

Closed-cycle refrigeration


Detector head = 40 lbs

Dark at 30V and  $-60^{\circ}\text{C}$  is .073 ADU/msec



2x Mode Photon Transfer Curve

- Read Noise = 2 ADU
- K constant =  $24 \text{ e-h/ADU} = 87 \text{ eV/ADU}$
- Full Well = 4000 ADU



8x Mode Photon Transfer Curve

- Read Noise = 7 ADU
- K constant =  $6 \text{ e-h/ADU} = 22 \text{ eV/ADU}$
- Full Well = 3500 ADU

# APS 480 x 480 Fast CCD Story - Future Tasks and Goals

- Convert GUI Software to EPICS Area Detector Software
  - Will allow spec to modify parameters through EPICS Channel Access
- Beam time at 8id on November 11, 2009
- Build more 480 x 480 Fast CCD Detector (1 at APS, 2 at ALS)
  - Build 4 more covers (1 for APS, 3 for ALS)
  - Test 3 Power Supplies (1 for APS, 2 for ALS)
  - Test out 3 Sets of Back End Electronics (1 set for APS, 2 sets for ALS)
- Maximize Frame Rate
  - Try fCRICs with faster ADC conversion time
    - Currently at 1.6usec per conversion which gives a 5.1msec readout
    - Would like to run fCRIC 1.0 usec per conversion or faster
  - Move system to Coreco X64 Xcelera-CL PX4 Frame Grabber
    - Uses XP 64 bit
    - Uses PCI Express x4 or x8 slot
    - Will allow us to increase frames per second for longer burst
  - Add faster Hard Drivers (Raid Array..) to maximize continuous frames per second
  - Try on-board background subtraction
  - Try on-board threshold data reduction



# Future of Fast CCD Detectors at ANL

- Building a 1k x 1k Frame Transfer Fast CCD
  - Using ARRA Funds
  - 2.5 Years to complete
  - Dedicated Wafer Run at Dalsa using LBNL's CCD Design
  - Using LBNL's fCRIC Chips
  - Scientist Requirements
    - Frame transfer
    - Smaller and lighter detector head, < 8 kg
    - Size of detector head < 200mm x 200mm x 200mm
    - Digitization  $\geq 1\text{MHz}$
    - On-board threshold and data compression



# The End

# Extra Slides




# The APS Detector Pool - Recent Additions

- Pilatus 100K digital Si pixel array detector developed at the Swiss Light Source/PSI (<http://pilatus.web.psi.ch/pilatus.htm>)
- a-Si flat panel detector (GE Healthcare; Perkin Elmer)
- 4-element Si drift detector array (SII NanoTechnology)




Pilatus 100K digital Si pixel array detector (SLS/PSI; XSD-BTS contact: A. Miceli)



a-Si flat panel detector (GE Healthcare; XSD-BTS contact: J. Lee)







4-element Si drift detector array (SII NanoTechnology; XSD-BTS contact: A. Miceli)




Same Photon transfer, but use image with same exposure as photon transfer images for average offset value.




# Dark Current Data at different Currents and Temperatures in 2x Mode



# Detector Improvements - Right to Left Intensity variation




# Detector Improvements - New 2x Configuration



Dark image with new 2x configuration

Same image  
-Over Scan Removed  
-B&C adjusted



Vertical cut of image  
top to bottom



# APS 480 x 480 Fast CCD Story - Science Data

## Scientific Motivation for Eye-Lens Protein Studies

- Understanding the glass transition in eye-lens protein mixtures will shed light on related physiological issues : in particular the stiffening of the eye lens seen in presbyopia, a loss of ability to focus the eye on nearby objects.
- Dynamic light scattering using lasers has measured a slowing down of long wavelength dynamics in mixtures of eye-lens proteins. However, DLS cannot see the protein-protein nearest neighbor peak, which should be the most sensitive to the onset of a glass transition.
- X-ray Photon Correlation Spectroscopy (XPCS) can be used to measure dynamics near this peak in concentrated suspensions. However, XPCS measurements on protein suspensions are extremely challenging due to fast dynamics and x-ray damage.
- The optimal camera for such experiments needs to work at high energy (to avoid x-ray damage to the proteins), have high efficiency (to make up for low signal) and have very fast readout (to maximize the range of time scales that can be

# APS 480 x 480 Fast CCD Story - Science Data

## Photon Droplet Binning - 11.5 keV

