

# Beamline 14-BM / BioCARS-CAT

**Scientific focus:** Structural biology

**Scientific programs:** Large unit cell (virus) crystallography, small unit cell (protein) crystallography, MAD phasing, time-resolved crystallography, Laue diffraction, and study of microcrystals

## Optics & Optical Performance

### 14-BM-C

- CARS-design conical focusing mirror  
25.731 m from source (on orbit along 3 mrad line)  
Si substrate, Rh coating  
water cooling at midplane  
4.1 mrad design angle  
horizontal focus: focusing monochromator  
vertical focus: bender
- CARS-design bent Ge(111) monochromator  
46.6 m from source (along 3 mrad centerline)  
Ge(111) single crystal  
bent triangle horizontal focusing
- misc. slits, collimators, filters, diagnostics etc.
- energy: capable 7.0–15.0 keV
- energy: operational 12.4 keV (1.0 A) fixed
- adjustable focal spot 210  $\mu\text{m}$  hor. x 350  $\mu\text{m}$  vert.  
(nominal)

### 14-BM-D

- 6.5–18.5 keV energy range, rapid tuning
- focal spot size 0.4 mm hor. x 0.6 mm vert. (nominal)  
adjustable to 150  $\mu\text{m}$  hor. x 150  $\mu\text{m}$  vert.
- CARS-design Si(111) double-bounce monochromator  
23.860 m from source (on orbit along 5 mrad line)  
6.5–18.5 energy range Si(111) crystal sets  
 $10^{-4}$  energy resolution ( $\Delta E/E$ ) at 10 keV  
38 mm nominal offset (fixed-exit, up bounce)  
water cooling
- CARS-design bent cylindrical focusing mirror  
Si substrate, Rh coating  
water cooling at midplane  
25.731 m from source (var. ht. along 5 mrad line)  
4.1 mrad design angle  
horizontal focus: sagittal cylindrical figure  
vertical focus: bender mechanism
- misc. slits, collimators, filters, diagnostics etc.

## Experiment Stations

### 14-BM-A

- white beam first optics enclosure
- 8.6 m x 1.8 m x 2.8 m (L x W x H)

### 14-BM-B

- pink beam optics enclosure
- 7.6 m x 1.2 m x 2.8 m (L x W x H)

### 14-BM-C

- monochromatic beam station
- 5.2 m x 3.7 m x 2.8 m (L x W x H)
- virus crystallography

### 14 BM-D

- pink or monochromatic beam station
- 6.4 m x 2.2 m x 2.8 m (L x W x H)
- MAD phasing
- time-resolved Laue (capable)
- protein crystallography

## Detectors

- 60° kappa diffractometer (all stations)
- single-axis diffractometer (available)
- ADSC Q1, ADSC Q4, MAR345 and off-line  
image plate detectors
- solid state, NaI scintillation, and Ge detectors

# Beamline 14-BM / BioCARS-CAT

**Scientific focus:** Structural biology

**Scientific programs:** Large unit cell (virus) crystallography, small unit cell (protein) crystallography, MAD phasing, time-resolved crystallography, Laue diffraction, and study of microcrystals

## Beamline Controls and Data Acquisition

- beamline and experimental control:
  - Dell Precision 620/PIII 933
  - Red Hat Linux 7.1
  - EPICS via VME
  - ADSC or Mar control software
- analysis: SGI 02 IRIX 6.5 running HKL/Denzo, DPS/MOSFLM, CCP4, CNS, LaueView, O, Predict, Resolve, Shelx, Strategy, Solve, XtalView, etc.

## Beamline Support Equipment/Facilities

- cryo-coolers: Oxford CryoStream, Oxford CryoJet, MSC, and CARS LN<sub>2</sub>/LHe<sub>2</sub> Cooler
- collimators, filters, slits, beam stop, CCD alignment cameras
- beam position monitors (1  $\mu\text{m}$  resolution)
- beam flux monitors
- BL3 facility, sample prep areas, cold room
- biochemistry equipment (pH meters, incubator, centrifuge, pipettes, glassware, lab refrigerator, etc.)

## Bending Magnet Source Characteristics (nominal)

source	APS bending magnet
critical energy	19.51 keV
on-axis peak brilliance at 16.3 keV	$2.9 \times 10^{15}$ ph/sec/mrad <sup>2</sup> /mm <sup>2</sup> /0.1%bw
on-axis peak angular flux at 16.3 keV	$9.6 \times 10^{13}$ ph/sec/mrad <sup>2</sup> /0.1%bw
on-axis peak horizontal angular flux at 5.6 keV	$1.6 \times 10^{13}$ ph/sec/mradh/0.1%bw
source size at critical energy $\sum_x$ $\sum_y$	145 $\mu\text{m}$ 36 $\mu\text{m}$
source divergence at critical energy $\sum_{x'}$ $\sum_{y'}$	6 mrad 47 $\mu\text{rad}$