areaDetector: A new module for EPICS area detector support

Mark Rivers

GeoSoilEnviroCARS, Advanced Photon Source

University of Chicago

areaDetector Talk Outline

- Brief overview of areaDetector module
- Highlight changes since last TWG presentation (June 19, 2008)
 - Releases
 - 1-0 (11-Apr-2008)
 - 1-1 (10-May-2008)
 - 1-2 (24-Oct-2008)
 - 1-3 (24-Nov-2008)
 - 1-4 (30-Jan-2009)
 - 1-5 (23-August-2009)
- Demo with Prosilica camera

areaDetector - Motivation


- 2-D detectors are essential components of synchrotron beamlines
 - Sample viewing cameras, x-ray diffraction and scattering detectors, xray imaging, optical spectroscopy, etc.
- Need to control the detectors from EPICS (useful even on non-EPICS beamlines, since other control systems like SPEC etc. can talk to EPICS)
- Previously several packages available, each typically restricted to a small set of detectors (Flea, Pilatus, marCCD, etc.)
- Clear advantages to an architecture that can be used on any detector, re-using many software components
- Providing EPICS control allows any higher-level client to control the detector and access the data (SPEC, medm, IDL programs, Python scripts, etc)

areaDetector - Goals

- Drivers for many detectors popular at synchrotron beamlines
 - Handle detectors ranging from >500 frames/second to <1 frame/second
- Basic parameters for all detectors
 - E.g. exposure time, start acquisition, etc.
 - Allows generic clients to be used for many applications
- Easy to implement new detector
 - Single device-driver C++ file to write. EPICS independent.
- Easy to implement detector-specific features
 - Driver understands additional parameters beyond those in the basic set
- EPICS-independent at lower layers.
- Middle-level plug-ins to add capability like regions-of-interest calculation, file saving, etc.
 - Device independent, work with all drivers
 - Below the EPICS layer for highest performance

areaDetector – Data structures

- NDArray
 - N-Dimensional array.
 - Everything is done in N-dimensions (up to 10), rather than 2. This is needed even for 2-D detectors to support color.
 - This is what plug-ins callbacks receive from device drivers.
- NDAttribute (NEW)
 - Each NDArray has a list of associated attributes (metadata) that travel with the array through the processing pileline. Attributes can come from driver parameters or any EPICS PV; e.g. can store motor positions, temperature, ring current, etc. with each frame.
- NDArrayPool
 - Allocates NDArray objects from a freelist
 - Plugins access in readonly mode, increment reference count
 - Eliminates need to copy data when sending it to callbacks.



Details: NDArray

```
/** Enumeration of NDArray data types */
typedef enum
{
   NDInt8, /**< Signed 8-bit integer */
   NDUInt8, /**< Unsigned 8-bit integer */
   NDInt16, /**< Signed 16-bit integer */
   NDUInt16, /**< Unsigned 16-bit integer */
   NDInt32, /**< Signed 32-bit integer */
   NDUInt32, /**< Unsigned 32-bit integer */
   NDFloat32, /**< 32-bit float */
   NDFloat64 /**< 64-bit float */
} NDDataType_t;
/** Structure defining a dimension of an NDArray */
typedef struct NDDimension {
   int size;
   int offset;
   int binning;
   int reverse;
} NDDimension t;
```

Details: NDArray

```
class NDArray {
•••
public:
    int
                  uniqueId;
    double
                  timeStamp;
    int
                  ndims;
    NDDimension_t dims[ND_ARRAY_MAX_DIMS];
    NDDataType_t dataType;
                  dataSize;
    int
    void
                  *pData;
    NDAttributeList *pAttributeList;
};
```

Details: NDAttribute

```
class NDAttribute {
public:
virtual int getValue(NDAttrDataType t dataType, void *pValue, size t
dataSize=0);
   virtual int setValue(NDAttrDataType t dataType, void *pValue);
   virtual int updateValue();
    ...
    char *pName;
    char *pDescription;
    char *pSource;
   NDAttrSource t sourceType;
   NDAttrDataType t dataType;
/** NDAttributeList class; this is a linked list of attributes. */
class NDAttributeList {
/** Attribute that gets its value from an EPICS PV. */
class PVAttribute : public NDAttribute {
/** Attribute that gets its value from an asynNDArrayDriver driver
parameter.
class paramAttribute : public NDAttribute {
```

Driver and PluginAttributes Defined via User-Created XML File

```
<?xml version="1.0" standalone="no" ?>
<Attributes>
  <Attribute name="AcquireTime"</pre>
             type="EPICS PV"
              source="13PS1:cam1:AcquireTime"
             dbrtype="DBR NATIVE"
             description="Camera acquire time" />
  <Attribute name="CameraModel"</pre>
             type="PARAM"
              source="MODEL"
             datatype="STRING"
             description="CameraModel" />
  <Attribute name="FramesDropped"</pre>
             type="PARAM"
              source="PS_FRAMES DROPPED"
             datatype="INT"
             description="FramesDropped" />
```

</Attributes>

Color Support (NEW)

- 4 colors models are supported
 - RGB1 (pixel interleave, e.g. RGB for each pixel are adjacent)
 - RGB2 (row interleave, e.g. R for row 1, then G, then B)
 - RGB3 (plane interleave, e.g. R for image, then G then B)
 - Bayer (color mask commonly used on color cameras)
- ColorConvert plugin converts between these color models
- Can reduce bandwidth 3X to read Bayer directly from camera (rather than RGB) and do color conversion in areaDetector

– Tradeoff bandwidth for CPU usage

• Can stream Bayer to file and do color conversion later, reduce file size 3X.

Detector drivers

- ADDriver
 - Base C++ class from which detector drivers derive. Handles details of EPICS interfaces, and other common functions.
- Simulation driver
 - Produces calculated images up to very high rates. Implements nearly all basic parameters, including color. Useful as a model for real detector drivers, and to test plugins and clients.
- Prosilica driver
 - Gigabit Ethernet cameras, mono and color
 - High resolution, high speed, e.g. 1360x1024 at 30 frames/second = 40MB/second.
- Firewire (IEEE-1396 DCAM) (NEW)
 - Vendor-independent Firewire camera drivers for Linux and Windows
- Roper driver (NEW)
 - Princeton Instruments and Photometrics cameras controlled via WinView

Detector drivers (continued)

- PVCAM driver (NEW)
 - Princeton Instruments and Photometrics cameras controlled via PVCAM library
- Pilatus driver (NEW)
 - Pilatus pixel-array detectors.
- marCCD driver (NEW)
 - Rayonix (MAR-USA) CCD x-ray detectors
- ADSC driver (NEW)
 - ADSC CCD detectors
- mar345 driver (NEW)
 - marResearch mar345 online image plate
- Perkin-Elmer driver (NEW)
 - Perkin-Elmer amorphous silicon detectors

ADBase.adl – Generic control screen

- Works with any detector
- Normally write custom control for each detector type to hide unimplemented features and expose driver-specific features

| Area Detector Control - 13SIM1:cam1:ShutterShutter mode Mome IShutter mode Mome IMonetedConnectedConnectedConnectedConnectedConnectedConnectedConnectedCollectKaposure time D.000 Close D.000PICS shutter setup DCollectExposure time D.000 0.010Acquire period D.000 0.010Acquire period D.000 0.010Acquire period D.000 0.010Acquire period D.000 0.010Mome IMome IImages complete 703# Exp./image IImage mode ContinuousInternalImage size 640 480Region size 640 480Region size 640 480Image size 640 | 🔀 ADBase.adl | | | | | | |
|--|--|--|--|--|--|--|--|
| asyn port SIM1 EPICS name 13SIM1:can1: Manufacturer Simulated detector Model Basic simulator Connected Connection Connect DisconnectShutter mode Mone Iss Status: Det. Closed EPICS Closed Open/Close Open Close Delay: Open 0.000 Close 0.000 EPICS shutter setup Iss Delay: Open 0.000 Close 0.000 # Images 10 10 # Image souplete 703 # Exp./image 1 1 Image mode Continuous Trigger mode Internal Internal Continuous Trigger mode Internal InternalRegion start 0000Region start 0000Region start 0000Region size 100480 Mo10001.00010001.00010001.000110001.000110001.000110001.0001110001.000 | Area Detector Control - 13SIM1:cam1: | | | | | | |
| EPICS name 13SIM1:cam1:Manufacturer Simulated detector Model Basic simulator ConnectedStatus: Det. Closed Open/Close Open CloseModel Basic simulator Connection ConnectDisconnectMore DDisconnectMore DDisconnectReadoutYSensor size 640480 1Binning II00Region start DD640480 480Region size 540480 480Region size 540480 480Region size 540480 480Reverse NoNoNoNoNoNoNoNoNoNoNoNoReverse NoNoNoNoData type Untes IUintes I | Setup | Shutter | | | | | |
| Manufacturer Simulated detector Model Basic simulator ConnectedOpen/Close Open Close Delay: Open 0.000ConnectedDelay: Open 0.000Close 0.000Connect DisconnectDisconnectCollectMore DCollectExposure time 0.0100.010ReadoutYYSensor size 640480# Images complete703IIIImage modeContinuousRegion start 000CollectorNoNoNoNoRegion size 640480Image counter 1Image size 640480Image counter 1NoNoNoNoNoNoNoNoImage rate 67.0Array callbacksEnableFileDriver file I/0Data type UnvsUInts | asyn port SIM1 | Shutter mode None | | | | | |
| Model Basic simulator ConnectedDelay: Open 0.000Close 0.000Connection Connect DisconnectEPICS shutter setup PMore PCollectMore PCollectMore PCollectMore PCollectMore PMore PReadoutXXYSensor size 64048011Binning 1100Region start 00640480Region size 540480Region size 540480Region size 540480Region size 540480Image size 640480Image size 640480Image size 640480Image size 640480Image size (bytes)307200Gain 10001.000Data type Unes 101Ulint8 | EPICS name 13SIM1:cam1: | Status: Det. Closed EPICS Closed | | | | | |
| Connected ConnectionEPICS shutter setupMoreDisconnectMoreCollectMoreExposure timeReadoutYXYSensor size64011Binning111Binning100Region start0640480Region size640800NoNoNoNoNoNoNoNoNoImage size640480Image size64011000Juage size100011000Data typeUints | Manufacturer Simulated detector | | | | | | |
| Lifts shutter setupConnectionConnectDisconnectMoreDisconnectCollectMoreDisconnectCollectReadoutXYSensor size64048011Image scomplete703BinningImageImage00Image modeContinuous00Image modeContinuousCollectingCollecting00Image modeContinuous11Image modeContinuous00Image modeContinuous00Image modeContinuous11Image modeContinuous00Image modeContinuous11Image modeContinuous11Image modeContinuous11Image modeContinuous11Image modeContinuous11Image modeContinuous11Image modeContinuous11Image modeContinuous11Image counter70311Image rate 67.0Image rate 67.011.0001.000File1Driver file I/0Di | | Delay: Open 0.000 Close 0.000 | | | | | |
| MoreCollectMoreCollectReadoutExposure timeXYSensor size64048011Binning100Region startP640480Region size640640480Region size640110000Region size640640480Region size6401111111111000000111 <t< th=""><th></th><th colspan="3">EPICS shutter setup 📃 🖻</th></t<> | | EPICS shutter setup 📃 🖻 | | | | | |
| MoreMo | Connection Connect Disconnect | Collect | | | | | |
| ReadoutAcquire period 0.0000.000ReadoutYImages 1010XYImages complete703Sensor size640480Image modeContinuousImageImageImageImageImage000ImageImage000ImageImage000ImageImage000ImageImage000ImageImage000ImageImage000ImageImage000ImageImage000ImageImage000ImageImage000ImageImage000ImageImage001.000ImageImage001.000ImageImage01.0001.000ImageImage001.000ImageImage001.000ImageImage001.000ImageImage001.000ImageImage001.000ImageImage001.000ImageImage001.000ImageImage001.000ImageImage001.000ImageImage0 </th <th>Mana D</th> <th>Exposure time 0.010 0.010</th> | Mana D | Exposure time 0.010 0.010 | | | | | |
| Readout# Images complete703XYSensor size640480111Binning11000Region start00640480Region size640640480Region size640MoNoNoNoReverseNoImage size640480Image counterImage size640Junage size640Mo1.000Image size640Gain1.000Junce1.000Data typeUInts | Hore 2 | Acquire period 0.000 0.000 | | | | | |
| XYSensor size64048011Binning111Binning100Region start000Region size640640480Region size640MoNoNoNoReverseNoImage size640480Image counterImage size64011 <t< th=""><th>Readout</th><th></th></t<> | Readout | | | | | | |
| Sensor size 640 480 1 1 Binning 1 0 0 Region start 0 640 480 Region size 640 640 480 Region size 640 No No No No Image size 640 Gain 1.000 Data type UInt8 | | | | | | | |
| 1 1 Binning 1 0 0 Region start 0 640 480 Region size 640 640 480 Region size 640 0 0 640 480 Region size 640 0 0 | | | | | | | |
| Binning I I I I I I I I I I I I I I I I I I I | | and a second | | | | | |
| Region start Image size 640 480 Region size 640 480 Detector state Image 0.000 No No No Image counter 703 Image size 640 480 Image rate 67.0 Image size (bytes) 307200 Array callbacks Enable Data type Imt8 UInt8 Driver file I/0 Image 1/0 | Binning 1 | Irigger mode Internal Internal | | | | | |
| 640 480 Region size 540 No No No No Image size 640 Image size 307200 Gain 1.000 Data type UInt8 | | Collecting | | | | | |
| Region size 640 480 Time remaining 0.000 No No No Image counter 0 703 Image size 640 480 Image rate 67.0 Array callbacks Enable Image size (bytes) 307200 Array callbacks Enable Enable Data type UInts UInts UInts Driver file I/0 D | and the second | | | | | | |
| No No Reverse No Image size 640 Image size 307200 Gain 1.000 Data type UInt8 | | | | | | | |
| Reverse No Image Image | | 5 | | | | | |
| Image Size Old Investigation Image Size (bytes) 307200 Gain 1.000 1.000 Data type UInts UInts | | | | | | | |
| Gain 1.000 1.000 Data type UInts UInts | | —————————————————————————————————————— | | | | | |
| Data type UInts UInt8 Driver file I/O | 0 5 | | | | | | |
| | | | | | | | |
| | | Uriver file 1/0 🔤 🖳 | | | | | |
| | COTOL IIIOGG LIOID | | | | | | |


Plugins

- Designed to perform real-time processing of data, running in the EPICS IOC (not over EPICS Channel Access)
- Receive NDArray data over callbacks from drivers or other plugins
- Plug-ins can execute in their own threads (non-blocking) or in callback thread (blocking)
 - If non-blocking then NDArray data is queued
 - Can drop images if queue is full
 - If executing in callback thread, no queuing, but slows device driver
- Allows
 - Enabling/disabling
 - Throttling rate (no more than 0.5 seconds, etc)
 - Changing data source for NDArray callbacks to another driver or plugin
- Some plugins are also sources of NDArray callbacks, as well as consumers.
 - Allows creating a data processing pipeline running at very high speed, each in a different thread, and hence in multiple cores on modern CPUs.

Plugins (continued)

- NDPlugInStdArrays
 - Receives arrays (images) from device drivers, converts to standard arrays, e.g. waveform records.
 - This plugin is what EPICS channel access viewers normally talk to.
- NDPluginROI
 - Performs region-of-interest calculations
 - Select a subregion. Optionally bin, reverse in either direction, convert data type. Optionally highlight borders within other ROIs.
 - Optionally compute statistics on the region (min, max, mean, total, net, etc.)
 - Optionally compute histogram of the region, make available as waveform for client plotting.
- NDPluginColorConvert (NEW)
 - Convert from one color model to another (Bayer, RGB pixel, row or planar interleave)
- NDPluginMJPEG (NEW)
 - MJPEG server that allows viewing images in a Web browser.

StdArrays and ROI plugin displays



Plugins: NDPluginFile

- Saves NDArrays to disk
- 3 modes:
 - Single array per disk file
 - Capture N arrays in memory, write to disk either multiple files or as a single large file (for file formats that support this.)
 - Stream arrays to a single large disk file
- File formats currently supported
 - TIFF (NEW)
 - JPEG (with compression control) (NEW)
 - netCDF (popular self-describing binary format, supported by Unidata at UCAR)
 - NeXus (NEW) (standard file format for neutron and x-ray communities, based on HDF, which is another popular self-describing binary format, richer than netCDF).
- netCDF and HDF save NDAttributes as well as NDArray data to the file
- In addition to file saving plugins, many vendor libraries also support saving files (e.g. marCCD, mar345, Pilatus, etc.) and this is supported at the driver level.
- File saving plugin can be used instead of or in addition to vendor file saving
 - Can add additional metadata vendor does not support
 - Could write JPEGS for Web display every minute, etc.

NDPluginFile display

| X NDFileNetCDF.adl | | | |
|--------------------|---|-----------------|---|
| | | 13SIM1:ne | tCDF1: |
| asyn port | SIM1FileNetCDF | | /home/epics/scratch/ |
| Plugin type | NDFileNetCDF | File path | /home/epics/scratch/ |
| Array port | BIM1 SIM1 | | test_netCDF |
| Array address | 0 | File name | test_netCDF |
| Enable | | Next file # | 14 14 |
| Min. time | the second se | Auto increment | Yes Yes |
| Callbacks block | | | %s%s_%d.nc |
| Array counter | 0 1695 | Filename format | %s%s_%d.nc File format netCDF InetCDF |
| Array rate | 80.0 | Last filename | /home/epics/scratch/test_netCDF_13.nc |
| Dropped arrays | 1015 C | | Hriting Done |
| # dimensions | 2 | Save file | Save Read file Read Auto save No P |
| Array Size | 640 480 0 | Write mode | stream - Stream # Capture 1000 1000 255 |
| Data type | UInt8 | | Capturing |
| Color mode | Mono | Capture | |
| Bayer pattern | RGGB | More | <u>-</u> |
| Unique ID | and the second se | | |
| | 619372961, 329 | | |
| Attributes file | petCDFAttributes.xml | | |

Example: streaming 80 frames/second of 640x480 video to netCDF file, no dropped frames.


Viewers

- areaDetector allows generic viewers to be written that receive images as EPICS waveform records over Channel Access
- Current viewers include:
 - ImageJ plugin EPICS_AD_Display. (NEW) ImageJ is a very popular image analysis program, written in Java, derived from NIH Image.
 - IDL EPICS_AD_Display.
 - IDL areaVision (NEW) (Stephen Mudie from AS)
 - MJPEG server (NEW) allows image display in any Web browser

ImageJ Viewer

| 🛓 ImageJ | | \mathbf{X} |
|---|---|--------------|
| File Edit Image Process Analyze Plugins Window Help | | |
| 🗆 🔾 📿 🗹 🔶 🔨 A 🔍 🖑 🌽 🕲 Dev. Stk. 🖉 🔏 🖄 | 1 | >> |
| x=610, y=421, value=86,92,89 | | |


| Image J EPICS_AD_Viewe | r Plugin | | | | | | | |
|------------------------|---------------|---------------|----------|----------|------------------|------|-------|------|
| PVPrefix | NX | NY | NZ | Frames/s | Capture to Stack | | | |
| 13PS1:image1: | 3 | 640 | 480 | 35.0 | | Snap | Start | Stop |
| Status: | 4/8/09 9:36:4 | 7.814: New im | nages=70 | | | | | |



Performance Example with Pilatus driver

- SPEC used to collect 1000 points using trajectory scanning mode with the Newport XPS motor controller. Hardware trigger of Pilatus from XPS.
- Relative scan of the chi axis from -2 degrees to +2 degrees with 1000 points at .02 seconds/point
- Coordinated motion of the phi, kappa and omega axes.
- Theoretical time 20.0 second, actual time 20.8 seconds
- Includes time to save all 1000 images to disk (366 MB), Pilatus driver to read each file, correct bad pixels and flat field, compute ROIs, and post the ROIs and

1000 images to EPICS.



Conclusions

- Architecture works well, easily extended to new detector drivers, new plugins and new clients
- Base classes, asynPortDriver, asynNDArrayDriver, asynPluginDriver actually are generic, nothing "areaDetector" specific about them.
- They can be used to implement any N-dimension detector, e.g. the XIA xMAP (16 detectors x 2048 channels x 512 points in a scan line)
- Now installed on all APS detector pool 2D detectors.
- Can get source code and pre-built binaries (Linux, Windows, Cygwin) from our Web site:
 - http://cars.uchicago.edu/software/epics/areaDetector
- Can also get code on SourceForge
 - http://epics.svn.sourceforge.net/viewvc/epics/applications/trunk/areaDetector

Acknowledgments

- Brian Tieman (APS) Roper PVCAM driver
- John Hammonds (APS) Perkin-Elmer driver and NeXus file saving plugin
- Tim Madden (APS) initial version of ImageJ viewer
- Stephen Mudie (Australian Synchrotron) areaViewer IDLbased viewer
- Ulrik Pedersen and Tom Cobb (Diamond) MJPEG plugin, Linux Firewire driver
- Lewis Muir (U Chicago) ADSC CCD driver
- NSF-EAR and DOE-Geosciences for support of GSECARS where most of this work was done
- Thanks for your attention!!!